

# BDP

Code: 0BDP

Disposable semiautomatic guillotine needle  
for coaxial biopsy on soft tissues

A



Choice between  
two different  
sampling depths

B



Releasable external  
cannula with  
Luer Lock connection

C



Atraumatic  
sharpened tip

D



High strength  
steel

E




Introducer needle  
(Trocar tip)  
Pre-position the needle for  
biopsy under TAC


F



Sliding stopper  
Internal echo-marker  
Depth marks



CE  
0426



# BDP

Disposable semiautomatic guillotine needle  
for coaxial biopsy on soft tissues

Code: 0BDP



M-72 I-51 I rev. A del 03/09/2012

To order:

CODE	SIZE	PACKAGE	COLOR C.
0BDP 1610	16Gx10 cm	20	WHITE
0BDP 1615	16Gx15 cm	20	WHITE
0BDP 1620	16Gx20 cm	20	WHITE
0BDP 1810	18Gx10 cm	20	PINK
0BDP 1815	18Gx15 cm	20	PINK
0BDP 1820	18Gx20 cm	20	PINK
0BDP 2015	20Gx15 cm	20	YELLOW
0BDP 2020	20Gx20 cm	20	YELLOW



On request Biomedical is available to provide needles with different lengths.

## biomedical

Biomedical srl Via G.B.Lulli, 43 50144 Firenze Italy  
Tel +39 055 352678 +39 055 332822 Fax +39 055 331629  
www.biomedical-bio.net sales@biomedical-bio.net

Bone biopsy needles

Guillotine needles for biopsy on soft tissues

Cyto-aspiration needles

Different specialistic needles and sets

Loco-regional anaesthesia needles and sets

Thermic and absorbing sheets