

Actionfold

Hydrophobic

Hydrophobic Aspheric Foldable Intraocular Lens



Model : AF-125600-PBA

CE
2460



Model : AF-125600-PBA
Hydrophobic Aspheric Optics

Hydrophobic Aspheric Foldable IOL

Salient Features :

- | | |
|--|---|
| <ul style="list-style-type: none"> • Type : Single Piece Acrylic Foldable • Overall length : 12.5 mm • Optic Diameter : 6.0 mm • Optic : Biconvex • Angulation : 0° (degree) • A-Constant : 118.5 • Refractive index : 1.481 • Sterilization : ETO STERILE EO • Power Range : +8.0D to +30.0D • UV Effect : UV Absorbing Hydrophobic | <ul style="list-style-type: none"> • Hydrophobic Material - Natural yellow imported from USA • 360° Square edge to minimize chances of PCO • Super intelligent memory for shape retention on unfolding • Bio-compatible hydro-acrylic composition • 6.0 mm optic enables to perform genuinely small incision surgery of less than 3.00 mm • Smooth folding & unfolding action • YAG Compatibility • High refractive index • Tear resistant • Natural yellow protects retina without blocking needed blue light • UV Blocking & Violet Blue Filtering |
|--|---|



An EN ISO 13485:2016 Certified Company

Action Medical Marketing Pvt. Ltd.

GAT No 528, Pune-Nagar Highway,
Near State Bank Of India, Koregaon Bhima, Tal. Shirur.

Dist: Pune 412 216. **INDIA**

Tel : 0091-2137-253617 / 18 / 19 / 20

Fax : 0091-2137-253620

Email : action@actionmedihealth.com

www.actionmedihealth.com

EC REP

Obelis S.A.
Bd. General Wahis 53,
B -1030 Brussels,
Belgium
E-mail : mail@obelis.net

HYDROPHOBIC 2 STEP PRELOADED LENS



SALIENT FEATURES

- Also Available in 2 step preloaded Delivery Systems
- Simple, effective, and provides ease of operation .
- Concentrate on the surgery, not on the delivery.



STERILE EO