

Envision

PMMA Rigid Intraocular Lens



Optic Style Equiconvex								
Model	EN-120500	EN-120525	EN-125550	EN-125602	EN-130602	EN-130652	EN-135652	EN-AC125
Overall length	12.0 mm	12.0 mm	12.0 mm	12.5 mm	13.0 mm	13.0 mm	13.5 mm	12.5 mm
Optic diameter	5.0 mm	5.25 mm	5.5 mm	6.0 mm	6.0 mm	6.5 mm	6.5 mm	6.0 mm
Positioning Holes	Nil	Nil	Nil	2	2	2	2	Nil
Hole position	-----	-----	-----	120 - 6'0 Clock	120 - 6'0 Clock	120 - 6'0 Clock	120 - 6'0 Clock	-----
A Constant	118.2	118.2	118.0	118.2	118.2	118.5	118.5	115.1
AC Depth	4.6 mm	4.6 mm	4.8 mm	4.8 mm	4.8 mm	5.0 mm	5.0 mm	2.9 mm

- ▶ **Choicest raw material**
- ▶ **Unique polishing process**
- ▶ **State-of-the-art technology**
- ▶ **Uncompromising quality control**
- ▶ **Choose - biconvex optics**
- ▶ **Ethylene oxide sterilization**
- ▶ **Selective models are also available in Square Edge (SQ) Optics.**

- ▶ **Premium quality IOL**
- ▶ **Super smooth optic surface finish**
- ▶ **Optimal combination of haptic strength : haptic flexibility**
- ▶ **Manufactured for operational & functional criticality**
- ▶ **Uses imported Tyvek pouches for assured sterility**
- ▶ **Internationally-styled lens inner case**
- ▶ **Material : PMMA (PolyMethylMethaAcrylate)**
- ▶ **Power range : -5.0D to +35.0D.**

Action

An EN ISO 13485:2016 Certified Company

Action Medical Marketing Pvt. Ltd.


GAT No 528, Pune-Nagar Highway, Near State Bank Of India,
Koregaon Bhima, Tal. Shirur, Dist. Pune 412 216. **INDIA**

Tel : 0091-2137-253617 / 18 / 19 / 20 Fax : 0091-2137-253620

Email : action@actionmedhealth.com www.actionmedhealth.com



Obelis S.A.
Bd. General Wafis 53,
B-1030 Brussels,
Belgium
E-mail : mail@obelis.net



Ver. 2/3/Sep. - 2019

■ Manufacturing

■ Marketing

■ Exporting