


# Hydrophilex

CE  
2460

## Acrylic Hydrophilic Aspheric Foldable Intraocular Lens




Model HF-125600



Model HF-125600

## Premium Hydrophilic Aspheric Foldable IOL

### Technical Specifications of Hydrophilex

- Type : Single Piece Acrylic Foldable
- Overall length : 12.5 mm
- Optic Diameter : 6.0 mm
- Optic : Equiconvex
- Angulations : 0° (degree)
- A-Constant : 118.2
- AC Depth : 5.2 mm
- Refractive index : 1.4626 (Hydrated)
- UV Effect : UV Absorbing
- Sterilization : Steam Sterilized 
- Power Range : -5.0D to +30.0D
- Dialing Slots : 0

### Salient Features for both Model (HF 125600 & HF 125600SQ)

- Hydrophilic Poly Hydroxy Ethyl Metha Acrylate (HEMA)
- Square edge to minimize change of PCO
- Super intelligent memory for shape retention on unfolding
- Bio-compatible hydro-acrylic composition
- 6.0 mm optic enables to perform genuinely small incision surgery of less than 3.00 mm
- Material has inherent UV absorbing capabilities
- Smooth folding & unfolding action
- YAG Compatibility
- High refractive index
- Tear resistant

## Action<sup>+</sup>

An EN ISO 13485:2016 Certified Company

**Action Medical Marketing Pvt. Ltd.**

GAT No 528, Pune-Nagar Highway,  
Near State Bank Of India, Koregaon Bhima, Tal. Shirur.

Dist: Pune 412 216. **INDIA**

Tel : 0091-2137-253617 / 18 / 19 / 20

Fax : 0091-2137-253620

Email : [action@actionmedihealth.com](mailto:action@actionmedihealth.com)

[www.actionmedihealth.com](http://www.actionmedihealth.com)

EC REP

Obelis S.A.  
Bd. General Wahis 53,  
B -1030 Brussels,  
Belgium  
E-mail : [mail@obelis.net](mailto:mail@obelis.net)

### Model HF-125600SQ

- Material : Hydrophilic Poly HEMA, Flexible, Biocompatible
- Dialing slot in haptic sector
- Optic Style : Equi - convex
- Optic Diameter : 6.0 mm
- Overall length : 12.5 mm
- "A" Constant : 118.2
- AC Depth : 5.2mm
- Diopter : -5.0D to +30.0D



Model HF-125600SQ

