

AMECATH

The logo for AMECATH features the word "AMECATH" in a bold, blue, sans-serif font. Below the text is a thin, light blue curved line that starts under the 'A' and ends under the 'H', curving upwards at both ends.

INTRAVASCULAR CATHETERS

A thick, dark blue wavy horizontal band that spans the width of the page, positioned below the main title.

AMECO MEDICAL INDUSTRIES

INTRAVASCULAR

Catheters and Accessories

1. Central Venous Catheters	3
1.1. Adult Central Venous Catheter Kits	4
1.2. Adult Pressure Resistant Central Venous Catheter Kits	11
1.3. Pediatric Central Venous Catheter Kits	14
1.4. Umbilical Catheter Kits	25
2. Peripheral Inserted Central Catheters (PICC)	27
2.1. Peripheral inserted Catheter kit	28
2.2. Pressure Resistant Peripheral Inserted Catheter kit	31
3. Midline Venous Catheters	33
3.1. Catheter Kit for Seldinger Technique (using IV Cath Cannula)	34
3.2. Midline Catheter Kit for Seldinger Technique (using over needle peel away introducer) :	35
3.3. Midline Catheter Kit for Seldinger Technique (using over needle peel away introducer) :	36
4. Arterial Catheters	37
4.1. Arterial catheters For Seldinger Technique kit	38
4.2. Arterial catheters For Seldinger Technique kit with Disposable Cartridge	38

AMECATH

CENTRAL VENOUS CATHETER KITS



1. Adult Central Venous Catheter Kits

KIT: CxLC-xxII-KN

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
5	15	CSLC- 0515-KN				
5	17	CSLC- 0517-KN				
5	20	CSLC- 0520-KN				
6	15	CSLC- 0615-KN	CDLC- 0615-KN	CTLC- 0615-KN		
6	17	CSLC- 0617-KN	CDLC- 0617-KN	CTLC- 0617-KN		
6	20	CSLC- 0620-KN	CDLC- 0620-KN	CTLC- 0620-KN		
7	15	CSLC- 0715-KN	CDLC- 0715-KN	CTLC- 0715-KN	CQLC- 0715-KN	CPLC- 0715-KN
7	17	CSLC- 0717-KN	CDLC- 0717-KN	CTLC- 0717-KN	CQLC- 0717-KN	CPLC- 0717-KN
7	20	CSLC- 0720-KN	CDLC- 0720-KN	CTLC- 0720-KN	CQLC- 0720-KN	CPLC- 0720-KN
7	30	CSLC- 0730-KN	CDLC- 0730-KN	CTLC- 0730-KN	CQLC- 0730-KN	CPLC- 0730-KN
7.5	20			CTLC- 07.520-KN		
8	15	CSLC- 0815-KN	CDLC- 0815-KN	CTLC- 0815-KN	CQLC- 0815-KN	CPLC- 0815-KN
8	17	CSLC- 0817-KN	CDLC- 0817-KN	CTLC- 0817-KN	CQLC- 0817-KN	CPLC- 0817-KN
8	20	CSLC- 0820-KN	CDLC- 0820-KN	CTLC- 0820-KN	CQLC- 0820-KN	CPLC- 0820-KN
8	30	CSLC- 0830-KN	CDLC- 0830-KN	CTLC- 0830-KN	CQLC- 0830-KN	CPLC- 0830-KN
8.5	15				CQLC- 08.515-KN	CPLC- 08.515-KN
8.5	17				CQLC- 08.517-KN	CPLC- 08.517-KN
8.5	20			CTLC- 08.520-KN	CQLC- 08.520-KN	CPLC- 08.520-KN
8.5	30				CQLC- 08.530-KN	CPLC- 08.530-KN
9	15				CQLC- 0915-KN	CPLC- 0915-KN
9	17				CQLC- 0917-KN	CPLC- 0917-KN
9	20				CQLC- 0920-KN	CPLC- 0920-KN
9	30				CQLC- 0930-KN	CPLC- 0930-KN

KIT CONTENTS :

- Adult Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing



KIT: CxLC-xxII-KN-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
5	15	CSLC- 0515-KN-A				
5	17	CSLC- 0517-KN-A				
5	20	CSLC- 0520-KN-A				
6	15	CSLC- 0615-KN-A	CDLC- 0615-KN-A	CTLC- 0615-KN-A		
6	17	CSLC- 0617-KN-A	CDLC- 0617-KN-A	CTLC- 0617-KN-A		
6	20	CSLC- 0620-KN-A	CDLC- 0620-KN-A	CTLC- 0620-KN-A		
7	15	CSLC- 0715-KN-A	CDLC- 0715-KN-A	CTLC- 0715-KN-A	CQLC- 0715-KN-A	CPLC- 0715-KN-A
7	17	CSLC- 0717-KN-A	CDLC- 0717-KN-A	CTLC- 0717-KN-A	CQLC- 0717-KN-A	CPLC- 0717-KN-A
7	20	CSLC- 0720-KN-A	CDLC- 0720-KN-A	CTLC- 0720-KN-A	CQLC- 0720-KN-A	CPLC- 0720-KN-A
7	30	CSLC- 0730-KN-A	CDLC- 0730-KN-A	CTLC- 0730-KN-A	CQLC- 0730-KN-A	CPLC- 0730-KN-A
7.5	20			CTLC- 07.520-KN-A		
8	15	CSLC- 0815-KN-A	CDLC- 0815-KN-A	CTLC- 0815-KN-A	CQLC- 0815-KN-A	CPLC- 0815-KN-A
8	17	CSLC- 0817-KN-A	CDLC- 0817-KN-A	CTLC- 0817-KN-A	CQLC- 0817-KN-A	CPLC- 0817-KN-A
8	20	CSLC- 0820-KN-A	CDLC- 0820-KN-A	CTLC- 0820-KN-A	CQLC- 0820-KN-A	CPLC- 0820-KN-A
8	30	CSLC- 0830-KN-A	CDLC- 0830-KN-A	CTLC- 0830-KN-A	CQLC- 0830-KN-A	CPLC- 0830-KN-A
8.5	15				CQLC- 08.515-KN-A	CPLC- 08.515-KN-A
8.5	17				CQLC- 08.517-KN-A	CPLC- 08.517-KN-A
8.5	20			CTLC- 08.520-KN-A	CQLC- 08.520-KN-A	CPLC- 08.520-KN-A
8.5	30				CQLC- 08.530-KN-A	CPLC- 08.530-KN-A
9	15				CQLC- 0915-KN-A	CPLC- 0915-KN-A
9	17				CQLC- 0917-KN-A	CPLC- 0917-KN-A
9	20				CQLC- 0920-KN-A	CPLC- 0920-KN-A
9	30				CQLC- 0930-KN-A	CPLC- 0930-KN-A

KIT CONTENTS :


- Adult Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing

KIT: CxLC-xxII-KGSN

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
5	15	CSLC- 0515-KGSN				
5	17	CSLC- 0517-KGSN				
5	20	CSLC- 0520-KGSN				
6	15	CSLC- 0615-KGSN	CDLC- 0615-KGSN	CTLC- 0615-KGSN		
6	17	CSLC- 0617-KGSN	CDLC- 0617-KGSN	CTLC- 0617-KGSN		
6	20	CSLC- 0620-KGSN	CDLC- 0620-KGSN	CTLC- 0620-KGSN		
7	15	CSLC- 0715-KGSN	CDLC- 0715-KGSN	CTLC- 0715-KGSN	CQLC- 0715-KGSN	CPLC- 0715-KGSN
7	16			CTLC- 0716-KGSN		
7	17	CSLC- 0717-KGSN	CDLC- 0717-KGSN	CTLC- 0717-KGSN	CQLC- 0717-KGSN	CPLC- 0717-KGSN
7	20	CSLC- 0720-KGSN	CDLC- 0720-KGSN	CTLC- 0720-KGSN	CQLC- 0720-KGSN	CPLC- 0720-KGSN
7	30	CSLC- 0730-KGSN	CDLC- 0730-KGSN	CTLC- 0730-KGSN	CQLC- 0730-KGSN	CPLC- 0730-KGSN
7.5	20			CTLC- 07.520-KGSN		
8	15	CSLC- 0815-KGSN	CDLC- 0815-KGSN	CTLC- 0815-KGSN	CQLC- 0815-KGSN	CPLC- 0815-KGSN
8	17	CSLC- 0817-KGSN	CDLC- 0817-KGSN	CTLC- 0817-KGSN	CQLC- 0817-KGSN	CPLC- 0817-KGSN
8	20	CSLC- 0820-KGSN	CDLC- 0820-KGSN	CTLC- 0820-KGSN	CQLC- 0820-KGSN	CPLC- 0820-KGSN
8	30	CSLC- 0830-KGSN	CDLC- 0830-KGSN	CTLC- 0830-KGSN	CQLC- 0830-KGSN	CPLC- 0830-KGSN
8.5	15			CTLC- 08.520-KGSN	CQLC- 08.515-KGSN	CPLC- 08.515-KGSN
8.5	17				CQLC- 08.517-KGSN	CPLC- 08.517-KGSN
8.5	20				CQLC- 08.520-KGSN	CPLC- 08.520-KGSN
8.5	30				CQLC- 08.530-KGSN	CPLC- 08.530-KGSN
9	15				CQLC- 0915-KGSN	CPLC- 0915-KGSN
9	17				CQLC- 0917-KGSN	CPLC- 0917-KGSN
9	20				CQLC- 0920-KGSN	CPLC- 0920-KGSN
9	30				CQLC- 0930-KGSN	CPLC- 0930-KGSN

KIT CONTENTS :

- Adult Central Venous Catheter
- Introducer needle
- Guiding Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing




KIT: CxLC-xxII-KGSN-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
5	15	CSLC- 0515-KGSN-A				
5	17	CSLC- 0517-KGSN-A				
5	20	CSLC- 0520-KGSN-A				
6	15	CSLC- 0615-KGSN-A	CDLC- 0615-KGSN-A	CTLC- 0615-KGSN-A		
6	17	CSLC- 0617-KGSN-A	CDLC- 0617-KGSN-A	CTLC- 0617-KGSN-A		
6	20	CSLC- 0620-KGSN-A	CDLC- 0620-KGSN-A	CTLC- 0620-KGSN-A		
7	15	CSLC- 0715-KGSN-A	CDLC- 0715-KGSN-A	CTLC- 0715-KGSN-A	CQLC- 0715-KGSN-A	CPLC- 0715-KGSN-A
7	16			CTLC- 0716-KGSN-A		
7	17	CSLC- 0717-KGSN-A	CDLC- 0717-KGSN-A	CTLC- 0717-KGSN-A	CQLC- 0717-KGSN-A	CPLC- 0717-KGSN-A
7	20	CSLC- 0720-KGSN-A	CDLC- 0720-KGSN-A	CTLC- 0720-KGSN-A	CQLC- 0720-KGSN-A	CPLC- 0720-KGSN-A
7	30	CSLC- 0730-KGSN-A	CDLC- 0730-KGSN-A	CTLC- 0730-KGSN-A	CQLC- 0730-KGSN-A	CPLC- 0730-KGSN-A
7.5	20			CTLC- 07.520-KGSN-A		
8	15	CSLC- 0815-KGSN-A	CDLC- 0815-KGSN-A	CTLC- 0815-KGSN-A	CQLC- 0815-KGSN-A	CPLC- 0815-KGSN-A
8	17	CSLC- 0817-KGSN-A	CDLC- 0817-KGSN-A	CTLC- 0817-KGSN-A	CQLC- 0817-KGSN-A	CPLC- 0817-KGSN-A
8	20	CSLC- 0820-KGSN-A	CDLC- 0820-KGSN-A	CTLC- 0820-KGSN-A	CQLC- 0820-KGSN-A	CPLC- 0820-KGSN-A
8	30	CSLC- 0830-KGSN-A	CDLC- 0830-KGSN-A	CTLC- 0830-KGSN-A	CQLC- 0830-KGSN-A	CPLC- 0830-KGSN-A
8.5	15			CTLC- 08.520-KGSN-A	CQLC- 08.515-KGSN-A	CPLC- 08.515-KGSN-A
8.5	17				CQLC- 08.517-KGSN-A	CPLC- 08.517-KGSN-A
8.5	20				CQLC- 08.520-KGSN-A	CPLC- 08.520-KGSN-A
8.5	30				CQLC- 08.530-KGSN-A	CPLC- 08.530-KGSN-A
9	15				CQLC- 0915-KGSN-A	CPLC- 0915-KGSN-A
9	17				CQLC- 0917-KGSN-A	CPLC- 0917-KGSN-A
9	20				CQLC- 0920-KGSN-A	CPLC- 0920-KGSN-A
9	30				CQLC- 0930-KGSN-A	CPLC- 0930-KGSN-A

KIT CONTENTS :

- Adult Central Venous Catheter
- Introducer needle
- Guiding Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing



KIT: CxLC-xxII-KGYN

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
5	15	CSLC- 0515-KGYN				
5	17	CSLC- 0517-KGYN				
5	20	CSLC- 0520-KGYN				
6	15	CSLC- 0615-KGYN	CDLC- 0615-KGYN	CTLC- 0615-KGYN		
6	17	CSLC- 0617-KGYN	CDLC- 0617-KGYN	CTLC- 0617-KGYN		
6	20	CSLC- 0620-KGYN	CDLC- 0620-KGYN	CTLC- 0620-KGYN		
7	15	CSLC- 0715-KGYN	CDLC- 0715-KGYN	CTLC- 0715-KGYN	CQLC- 0715-KGYN	CPLC- 0715-KGYN
7	16					
7	17	CSLC- 0717-KGYN	CDLC- 0717-KGYN	CTLC- 0717-KGYN	CQLC- 0717-KGYN	CPLC- 0717-KGYN
7	20	CSLC- 0720-KGYN	CDLC- 0720-KGYN	CTLC- 0720-KGYN	CQLC- 0720-KGYN	CPLC- 0720-KGYN
7	30	CSLC- 0730-KGYN	CDLC- 0730-KGYN	CTLC- 0730-KGYN	CQLC- 0730-KGYN	CPLC- 0730-KGYN
7.5	20					
8	15	CSLC- 0815-KGYN	CDLC- 0815-KGYN	CTLC- 0815-KGYN	CQLC- 0815-KGYN	CPLC- 0815-KGYN
8	17	CSLC- 0817-KGYN	CDLC- 0817-KGYN	CTLC- 0817-KGYN	CQLC- 0817-KGYN	CPLC- 0817-KGYN
8	20	CSLC- 0820-KGYN	CDLC- 0820-KGYN	CTLC- 0820-KGYN	CQLC- 0820-KGYN	CPLC- 0820-KGYN
8	30	CSLC- 0830-KGYN	CDLC- 0830-KGYN	CTLC- 0830-KGYN	CQLC- 0830-KGYN	CPLC- 0830-KGYN
8.5	15				CQLC- 08.515-KGYN	CPLC- 08.515-KGYN
8.5	17				CQLC- 08.517-KGYN	CPLC- 08.517-KGYN
8.5	20				CQLC- 08.520-KGYN	CPLC- 08.520-KGYN
8.5	30				CQLC- 08.530-KGYN	CPLC- 08.530-KGYN
9	15				CQLC- 0915-KGYN	CPLC- 0915-KGYN
9	17				CQLC- 0917-KGYN	CPLC- 0917-KGYN
9	20				CQLC- 0920-KGYN	CPLC- 0920-KGYN
9	30				CQLC- 0930-KGYN	CPLC- 0930-KGYN

KIT CONTENTS :


- Adult Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Guiding Y connector
- Fixation Wing

KIT: CxLC-xxII-KGYN-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
5	15	CSLC- 0515-KGYN-A				
5	17	CSLC- 0517-KGYN-A				
5	20	CSLC- 0520-KGYN-A				
6	15	CSLC- 0615-KGYN-A	CDLC- 0615-KGYN-A	CTLC- 0615-KGYN-A		
6	17	CSLC- 0617-KGYN-A	CDLC- 0617-KGYN-A	CTLC- 0617-KGYN-A		
6	20	CSLC- 0620-KGYN-A	CDLC- 0620-KGYN-A	CTLC- 0620-KGYN-A		
7	15	CSLC- 0715-KGYN-A	CDLC- 0715-KGYN-A	CTLC- 0715-KGYN-A	CQLC- 0715-KGYN-A	CPLC- 0715-KGYN-A
7	16					
7	17	CSLC- 0717-KGYN-A	CDLC- 0717-KGYN-A	CTLC- 0717-KGYN-A	CQLC- 0717-KGYN-A	CPLC- 0717-KGYN-A
7	20	CSLC- 0720-KGYN-A	CDLC- 0720-KGYN-A	CTLC- 0720-KGYN-A	CQLC- 0720-KGYN-A	CPLC- 0720-KGYN-A
7	30	CSLC- 0730-KGYN-A	CDLC- 0730-KGYN-A	CTLC- 0730-KGYN-A	CQLC- 0730-KGYN-A	CPLC- 0730-KGYN-A
7.5	20					
8	15	CSLC- 0815-KGYN-A	CDLC- 0815-KGYN-A	CTLC- 0815-KGYN-A	CQLC- 0815-KGYN-A	CPLC- 0815-KGYN-A
8	17	CSLC- 0817-KGYN-A	CDLC- 0817-KGYN-A	CTLC- 0817-KGYN-A	CQLC- 0817-KGYN-A	CPLC- 0817-KGYN-A
8	20	CSLC- 0820-KGYN-A	CDLC- 0820-KGYN-A	CTLC- 0820-KGYN-A	CQLC- 0820-KGYN-A	CPLC- 0820-KGYN-A
8	30	CSLC- 0830-KGYN-A	CDLC- 0830-KGYN-A	CTLC- 0830-KGYN-A	CQLC- 0830-KGYN-A	CPLC- 0830-KGYN-A
8.5	15				CQLC- 08.515-KGYN-A	CPLC- 08.515-KGYN-A
8.5	17				CQLC- 08.517-KGYN-A	CPLC- 08.517-KGYN-A
8.5	20				CQLC- 08.520-KGYN-A	CPLC- 08.520-KGYN-A
8.5	30				CQLC- 08.530-KGYN-A	CPLC- 08.530-KGYN-A
9	15				CQLC- 0915-KGYN-A	CPLC- 0915-KGYN-A
9	17				CQLC- 0917-KGYN-A	CPLC- 0917-KGYN-A
9	20				CQLC- 0920-KGYN-A	CPLC- 0920-KGYN-A
9	30				CQLC- 0930-KGYN-A	CPLC- 0930-KGYN-A

KIT CONTENTS :

- Adult Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Guiding Y connector
- Fixation Wing



KIT: CxLC-xxII-KGSNO

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
5	15	CSLC- 0515-KGSNO		
5	17	CSLC- 0517-KGSNO		
5	20	CSLC- 0520-KGSNO		
6	15	CSLC- 0615-KGSNO	CDLC- 0615-KGSNO	CTLC- 0615-KGSNO
6	17	CSLC- 0617-KGSNO	CDLC- 0617-KGSNO	CTLC- 0617-KGSNO
6	20	CSLC- 0620-KGSNO	CDLC- 0620-KGSNO	CTLC- 0620-KGSNO
7	15	CSLC- 0715-KGSNO	CDLC- 0715-KGSNO	CTLC- 0715-KGSNO
7	16			
7	17	CSLC- 0717-KGSNO	CDLC- 0717-KGSNO	CTLC- 0717-KGSNO
7	20	CSLC- 0720-KGSNO	CDLC- 0720-KGSNO	CTLC- 0720-KGSNO
7	30	CSLC- 0730-KGSNO	CDLC- 0730-KGSNO	CTLC- 0730-KGSNO
7	60		CDLC-0760-KGSNO	CTLC- 0760-KGSNO
7.5	20			
8	15	CSLC- 0815-KGSNO	CDLC- 0815-KGSNO	CTLC- 0815-KGSNO
8	17	CSLC- 0817-KGSNO	CDLC- 0817-KGSNO	CTLC- 0817-KGSNO
8	20	CSLC- 0820-KGSNO	CDLC- 0820-KGSNO	CTLC- 0820-KGSNO
8	30	CSLC- 0830-KGSNO	CDLC- 0830-KGSNO	CTLC- 0830-KGSNO
8.5	15			
8.5	17			
8.5	20			
8.5	30			
9	15			
9	17			
9	20			

KIT CONTENTS :

- Adult Central Venous Catheter
- Over the needle catheter
- Introducer needle
- Guiding Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection cap
- Transparent catheter fixation adhesive

KIT: CxLC-xxII-KGSNO-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
5	15	CSLC- 0515-KGSNO-A		
5	17	CSLC- 0517-KGSNO-A		
5	20	CSLC- 0520-KGSNO-A		
6	15	CSLC- 0615-KGSNO-A	CDLC- 0615-KGSNO-A	CTLC- 0615-KGSNO-A
6	17	CSLC- 0617-KGSNO-A	CDLC- 0617-KGSNO-A	CTLC- 0617-KGSNO-A
6	20	CSLC- 0620-KGSNO-A	CDLC- 0620-KGSNO-A	CTLC- 0620-KGSNO-A
7	15	CSLC- 0715-KGSNO-A	CDLC- 0715-KGSNO-A	CTLC- 0715-KGSNO-A
7	16			
7	17	CSLC- 0717-KGSNO-A	CDLC- 0717-KGSNO-A	CTLC- 0717-KGSNO-A
7	20	CSLC- 0720-KGSNO-A	CDLC- 0720-KGSNO-A	CTLC- 0720-KGSNO-A
7	30	CSLC- 0730-KGSNO-A	CDLC- 0730-KGSNO-A	CTLC- 0730-KGSNO-A
7	60		CDLC-0760-KGSNO-A	CTLC- 0760-KGSNO-A
7.5	20			
8	15	CSLC- 0815-KGSNO-A	CDLC- 0815-KGSNO-A	CTLC- 0815-KGSNO-A
8	17	CSLC- 0817-KGSNO-A	CDLC- 0817-KGSNO-A	CTLC- 0817-KGSNO-A
8	20	CSLC- 0820-KGSNO-A	CDLC- 0820-KGSNO-A	CTLC- 0820-KGSNO-A
8	30	CSLC- 0830-KGSNO-A	CDLC- 0830-KGSNO-A	CTLC- 0830-KGSNO-A
8.5	15			
8.5	17			
8.5	20			
8.5	30			
9	15			
9	17			
9	20			
9	30			

KIT CONTENTS :

- Adult Central Venous Catheter
- Over the needle catheter
- Introducer needle
- Guiding Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection cap
- Transparent catheter fixation adhesive

2. Adult Pressure Resistant Central Venous Catheter Kits

KIT:PRCxLC-xxII-KN

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
6	15	PRCSLC-0615-KN	PRCDLC-0615-KN	PRCTLTC-0615-KN		
6	20	PRCSLC-0620-KN	PRCDLC-0620-KN	PRCTLTC-0620-KN		
7	15	PRCSLC-0715-KN	PRCDLC-0715-KN	PRCTLTC-0715-KN		
7	20	PRCSLC-0720-KN	PRCDLC-0720-KN	PRCTLTC-0720-KN		
8	15	PRCSLC-0815-KN	PRCDLC-0815-KN	PRCTLTC-0815-KN	PRCQLC- 0815-KN	PRCPLC- 0815-KN
8	20	PRCSLC-0820-KN	PRCDLC-0820-KN	PRCTLTC-0820-KN	PRCQLC- 0820-KN	PRCPLC- 0820-KN
9	15				PRCQLC- 0915-KN	PRCPLC- 0915-KN
9	20				PRCQLC- 0920-KN	PRCPLC- 0920-KN

KIT CONTENTS :

- Adult Pressure Resistant Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing

KIT:PRCxLC-xxII-KN-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
6	15	PRCSLC-0615-KN-A	PRCDLC-0615-KN-A	PRCTLTC-0615-KN-A		
6	20	PRCSLC-0620-KN-A	PRCDLC-0620-KN-A	PRCTLTC-0620-KN-A		
7	15	PRCSLC-0715-KN-A	PRCDLC-0715-KN-A	PRCTLTC-0715-KN-A		
7	20	PRCSLC-0720-KN-A	PRCDLC-0720-KN-A	PRCTLTC-0720-KN-A		
8	15	PRCSLC-0815-KN-A	PRCDLC-0815-KN-A	PRCTLTC-0815-KN-A	PRCQLC- 0815-KN-A	PRCPLC- 0815-KN-A
8	20	PRCSLC-0820-KN-A	PRCDLC-0820-KN-A	PRCTLTC-0820-KN-A	PRCQLC- 0820-KN-A	PRCPLC- 0820-KN-A
9	15				PRCQLC- 0915-KN-A	PRCPLC- 0915-KN-A
9	20				PRCQLC- 0920-KN-A	PRCPLC- 0920-KN-A

KIT CONTENTS :

- Adult Pressure Resistant Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing

KIT:PRCxLC-xxII-KGSN

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
6	15	PRCSLC-0615-KGSN	PRCDLC- 0615-KGSN	PRCTL- 0615-KGSN		
6	20	PRCSLC-0620-KGSN	PRCDLC- 0620- KGSN	PRCTL- 0620- KGSN		
7	15	PRCSLC-0715-KGSN	PRCDLC- 0715- KGSN	PRCTL- 0715- KGSN		
7	20	PRCSLC-0720-KGSN	PRCDLC- 0720- KGSN	PRCTL- 0720- KGSN		
8	15	PRCSLC-0815-KGSN	PRCDLC- 0815- KGSN	PRCTL- 0815- KGSN	PRCQLC- 0815-KGSN	PRCPLC- 0815-KGSN
8	20	PRCSLC-0820-KGSN	PRCDLC- 0820- KGSN	PRCTL- 0820- KGSN	PRCQLC- 0820-KGSN	PRCPLC- 0820-KGSN
9	15				PRCQLC- 0915-KGSN	PRCPLC- 0915-KGSN
9	20				PRCQLC- 0920-KGSN	PRCPLC- 0920-KGSN

KIT CONTENTS :

- Adult Pressure Resistant Central Venous Catheter
- Introducer needle
- Guiding Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing

KIT:PRCxLC-xxII-KGSN-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
6	15	PRCSLC-0615-KGSN-A	PRCDLC- 0615-KGSN-A	PRCTL- 0615-KGSN-A		
6	20	PRCSLC-0620-KGSN-A	PRCDLC- 0620- KGSN-A	PRCTL- 0620- KGSN-A		
7	15	PRCSLC-0715-KGSN-A	PRCDLC- 0715- KGSN-A	PRCTL- 0715- KGSN-A		
7	20	PRCSLC-0720-KGSN-A	PRCDLC- 0720- KGSN-A	PRCTL- 0720- KGSN-A		
8	15	PRCSLC-0815-KGSN-A	PRCDLC- 0815- KGSN-A	PRCTL- 0815- KGSN-A	PRCQLC- 0815-KGSN-A	PRCPLC- 0815-KGSN-A
8	20	PRCSLC-0820-KGSN-A	PRCDLC- 0820- KGSN-A	PRCTL- 0820- KGSN-A	PRCQLC- 0820-KGSN-A	PRCPLC- 0820-KGSN-A
9	15				PRCQLC- 0915-KGSN-A	PRCPLC- 0915-KGSN-A
9	20				PRCQLC- 0920-KGSN-A	PRCPLC- 0920-KGSN-A

KIT CONTENTS :

- Adult Pressure Resistant Central Venous Catheter
- Introducer needle
- Guiding Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing



KIT:PRCxLC-xxII-KGYN

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
6	15	PRCSLC-0615-KGYN	PRCDLC- 0615-KGYN	PRCTLC- 0615-KGYN		
6	20	PRCSLC-0620-KGYN	PRCDLC- 0620- KGYN	PRCTLC- 0620- KGYN		
7	15	PRCSLC-0715-KGYN	PRCDLC- 0715- KGYN	PRCTLC- 0715- KGYN		
7	20	PRCSLC-0720-KGYN	PRCDLC- 0720- KGYN	PRCTLC- 0720- KGYN		
8	15	PRCSLC-0815-KGYN	PRCDLC- 0815- KGYN	PRCTLC- 0815- KGYN	PRCQLC- 0815-KGYN	PRCPLC- 0815-KGYN
8	20	PRCSLC-0820-KGYN		PRCTLC- 0820- KGYN	PRCQLC- 0820-KGYN	PRCPLC- 0820-KGYN
9	15				PRCQLC- 0915-KGYN	PRCPLC- 0915-KGYN
9	20				PRCQLC- 0920-KGYN	PRCPLC- 0920-KGYN

KIT CONTENTS :

- Adult Pressure Resistant Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Guiding Y connector
- Fixation Wing

KIT:PRCxLC-xxII-KGYN-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen	Penta lumen
6	15	PRCSLC-0615-KGYN-A	PRCDLC-0615-KGYN-A	PRCTLC- 0615-KGYN-A		
6	20	PRCSLC-0620-KGYN-A	PRCDLC-0620-KGYN-A	PRCTLC- 0620- KGYN-A		
7	15	PRCSLC-0715-KGYN-A	PRCDLC-0715-KGYN-A	PRCTLC- 0715- KGYN-A		
7	20	PRCSLC-0720-KGYN-A	PRCDLC-0720-KGYN-A	PRCTLC- 0720- KGYN-A		
8	15	PRCSLC-0815-KGYN-A	PRCDLC-0815-KGYN-A	PRCTLC- 0815- KGYN-A	PRCQLC- 0815-KGYN-A	PRCPLC- 0815-KGYN-A
8	20	PRCSLC-0820-KGYN-A		PRCTLC- 0820- KGYN-A	PRCQLC- 0820-KGYN-A	PRCPLC- 0820-KGYN-A
9	15				PRCQLC- 0915-KGYN-A	PRCPLC- 0915-KGYN-A
9	20				PRCQLC- 0920-KGYN-A	PRCPLC- 0920-KGYN-A

KIT CONTENTS :

- Adult Pressure Resistant Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Guiding Y connector
- Fixation Wing

3. Pediatric Central Venous Catheter Kits

KIT:CxLC-xxII-KP

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
2	5	CSLC- 0205-KP		
2	8	CSLC- 0208-KP		
3	4		CDLC- 0304-KP	
3	5	CSLC- 0305-KP	CDLC- 0305-KP	
3	8	CSLC- 0308-KP	CDLC- 0308-KP	
3	10	CSLC- 0310-KP	CDLC- 0310-KP	
3	13	CSLC- 0313-KP	CDLC- 0313-KP	
3	15	CSLC- 0315-KP		
3	17	CSLC- 0317-KP		
4	5	CSLC- 0405-KP	CDLC- 0405-KP	CTLC- 0405-KP
4	6	CSLC- 0406-KP	CDLC- 0406-KP	CTLC- 0406-KP
4	8	CSLC- 0408-KP	CDLC- 0408-KP	CTLC- 0408-KP
4	10	CSLC- 0410-KP	CDLC- 0410-KP	CTLC- 0410-KP
4	12	CSLC- 0412-KP		
4	13	CSLC- 0413-KP	CDLC- 0413-KP	CTLC- 0413-KP
4	15	CSLC- 0415-KP	CDLC- 0415-KP	
4	17	CSLC- 0417-KP		
4	20	CSLC- 0420-KP		
4	30			CTLC- 0430-KP
4.5	8		CDLC- 04.508-KP	
4.5	10		CDLC- 04.510-KP	



KIT:CxLC-xxII-KP

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen
5	5		CDLC- 0505-KP	CTLC- 0505-KP	
5	6	CSLC- 0506-KP	CDLC- 0506-KP		
5	8	CSLC- 0508-KP	CDLC- 0508-KP	CTLC- 0508-KP	
5	10	CSLC- 0510-KP	CDLC- 0510-KP	CTLC- 0510-KP	
5	13	CSLC- 0513-KP	CDLC- 0513-KP	CTLC- 0513-KP	
5	15	CSLC- 0515-KP	CDLC- 0515-KP	CTLC- 0515-KP	
5	17	CSLC- 0517-KP			
5	20	CSLC- 0520-KP		CTLC- 0520-KP	
5	30			CTLC- 0530-KP	
5.5	10		CDLC- 05.510-KP	CTLC- 05.510-KP	
5.5	15		CDLC- 05.515-KP	CTLC- 05.515-KP	
6	8				CQLC- 0608-KP
6	13				CQLC- 0613-KP
6	15				CQLC- 0615-KP
6	20				CQLC- 0620-KP

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection caps



KIT:CxLC-xxII-KP-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
2	5	CSLC- 0205-KP-A		
2	8	CSLC- 0208-KP-A		
3	4		CDLC- 0304-KP-A	
3	5	CSLC- 0305-KP-A	CDLC- 0305-KP-A	
3	8	CSLC- 0308-KP-A	CDLC- 0308-KP-A	
3	10	CSLC- 0310-KP-A	CDLC- 0310-KP-A	
3	13	CSLC- 0313-KP-A	CDLC- 0313-KP-A	
3	15	CSLC- 0315-KP-A		
3	17	CSLC- 0317-KP-A		
4	5	CSLC- 0405-KP-A	CDLC- 0405-KP-A	CTLC- 0405-KP-A
4	6	CSLC- 0406-KP-A	CDLC- 0406-KP-A	CTLC- 0406-KP-A
4	8	CSLC- 0408-KP-A	CDLC- 0408-KP-A	CTLC- 0408-KP-A
4	10	CSLC- 0410-KP-A	CDLC- 0410-KP-A	CTLC- 0410-KP-A
4	12	CSLC- 0412-KP-A		
4	13	CSLC- 0413-KP-A	CDLC- 0413-KP-A	CTLC- 0413-KP-A
4	15	CSLC- 0415-KP-A	CDLC- 0415-KP-A	
4	17	CSLC- 0417-KP-A		
4	20	CSLC- 0420-KP-A		
4	30			CTLC- 0430-KP-A
4.5	8		CDLC- 04.508-KP-A	
4.5	10		CDLC- 04.510-KP-A	



KIT:CxLC-xxII-KP

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen
5	5		CDLC- 0505-KP-A	CTLC- 0505-KP-A	
5	6	CSLC- 0506-KP-A	CDLC- 0506-KP-A		
5	8	CSLC- 0508-KP-A	CDLC- 0508-KP-A	CTLC- 0508-KP-A	
5	10	CSLC- 0510-KP-A	CDLC- 0510-KP-A	CTLC- 0510-KP-A	
5	13	CSLC- 0513-KP-A	CDLC- 0513-KP-A	CTLC- 0513-KP-A	
5	15	CSLC- 0515-KP-A	CDLC- 0515-KP-A	CTLC- 0515-KP-A	
5	17	CSLC- 0517-KP-A			
5	20	CSLC- 0520-KP-A		CTLC- 0520-KP-A	
5	30			CTLC- 0530-KP-A	
5.5	10		CDLC- 05.510-KP-A	CTLC- 05.510-KP-A	
5.5	15		CDLC- 05.515-KP-A	CTLC- 05.515-KP-A	
6	8				CQLC- 0608-KP-A
6	13				CQLC- 0613-KP-A
6	15				CQLC- 0615-KP-A
6	20				CQLC- 0620-KP-A

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection Caps

KIT:CxLC-xxII-KP*

Fr. Size	Length	Single lumen
2	5	CSLC- 0205-KP*
2	8	CSLC- 0208-KP*
2	15	CSLC- 0215-KP*
2	30	CSLC- 0230-KP*
3	5	CSLC- 0305-KP*
3	8	CSLC- 0308-KP*
3	13	CSLC- 0313-KP*
3	15	CSLC- 0315-KP*
3	30	CSLC- 0330-KP*

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Peel Away over the Needle
- Measuring tape
- T-connector With Stiffener Wire

KIT:CxLC-xxII-KP**

Fr. Size	Length	Single lumen	Dual lumen
3	5	CSLC- 0305-KP**	CDLC- 0305-KP**
3	8	CSLC- 0308-KP**	CDLC- 0308-KP**
3	13	CSLC- 0313-KP**	CDLC- 0313-KP**
3	15	CSLC- 0315-KP**	CDLC- 0315-KP**
3	30	CSLC- 0330-KP**	CDLC- 0330-KP**

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer Needles
- Peel away over Dilator
- Syringe
- Nitinol guide wire J end with Dispenser
- Measuring tape
- T-connector with Stiffener Wire

KIT:CxLC-xxII-KP**-A

Fr. Size	Length	Single lumen	Dual lumen
3	5	CSLC- 0305-KP**-A	CDLC- 0305-KP**-A
3	8	CSLC- 0308-KP**-A	CDLC- 0308-KP**-A
3	13	CSLC- 0313-KP**-A	CDLC- 0313-KP**-A
3	15	CSLC- 0315-KP**-A	CDLC- 0315-KP**-A
3	30	CSLC- 0330-KP**-A	CDLC- 0330-KP**-A

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer Needles
- Peel away over Dilator
- Syringe
- Nitinol guide wire J end with Advancer
- Measuring tape
- T-connector with Stiffener Wire

KIT:CxLC-xxII-KPO

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen
2	5	CSLC- 0205-KPO			
2	8	CSLC- 0208-KPO			
3	4		CDLC- 0304-KPO		
3	5	CSLC- 0305-KPO	CDLC- 0305-KPO		
3	8	CSLC- 0308-KPO	CDLC- 0308-KPO		
3	10	CSLC- 0310-KPO	CDLC- 0310-KPO		
3	13	CSLC- 0313-KPO	CDLC- 0313-KPO		
3	15	CSLC- 0315-KPO			
3	17	CSLC- 0317-KPO			
4	5	CSLC- 0405-KPO	CDLC- 0405-KPO	CTLC- 0405-KPO	
4	6	CSLC- 0406-KPO	CDLC- 0406-KPO	CTLC- 0406-KPO	
4	8	CSLC- 0408-KPO	CDLC- 0408-KPO	CTLC- 0408-KPO	
4	10	CSLC- 0410-KPO	CDLC- 0410-KPO	CTLC- 0410-KPO	
4	12	CSLC- 0412-KPO			
4	13	CSLC- 0413-KPO	CDLC- 0413-KPO	CTLC- 0413-KPO	
4	15	CSLC- 0415-KPO	CDLC- 0415-KPO		
4	17	CSLC- 0417-KPO			
4	20	CSLC- 0420-KPO			
4	30			CTLC- 0430-KPO	
4.5	8		CDLC- 04.508-KPO		
4.5	10		CDLC- 04.510-KPO		
5	5		CDLC- 0505-KPO	CTLC- 0505-KPO	
5	6	CSLC- 0506-KPO	CDLC- 0506-KPO		
5	8	CSLC- 0508-KPO	CDLC- 0508-KPO	CTLC- 0508-KPO	
5	10	CSLC- 0510-KPO	CDLC- 0510-KPO	CTLC- 0510-KPO	
5	13	CSLC- 0513-KPO	CDLC- 0513-KPO	CTLC- 0513-KPO	
5	15	CSLC- 0515-KPO	CDLC- 0515-KPO	CTLC- 0515-KPO	
5	17	CSLC- 0517-KPO			
5	20	CSLC- 0520-KPO		CTLC- 0520-KPO	
5	30			CTLC- 0530-KPO	
5.5	10		CDLC- 05.510-KPO	CTLC- 05.510-KPO	
5.5	15		CDLC- 05.515-KPO	CTLC- 05.515-KPO	
6	8				CQLC- 0608-KPO
6	13				CQLC- 0613-KPO
6	15				CQLC- 0615-KPO
6	20				CQLC- 0620-KPO

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection caps

KIT:CxLC-xxII-KPO-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen
2	5	CSLC- 0205-KPO-A			
2	8	CSLC- 0208-KPO-A			
3	4		CDLC- 0304-KPO-A		
3	5	CSLC- 0305-KPO-A	CDLC- 0305-KPO-A		
3	8	CSLC- 0308-KPO-A	CDLC- 0308-KPO-A		
3	10	CSLC- 0310-KPO-A	CDLC- 0310-KPO-A		
3	13	CSLC- 0313-KPO-A	CDLC- 0313-KPO-A		
3	15	CSLC- 0315-KPO-A			
3	17	CSLC- 0317-KPO-A			
4	5	CSLC- 0405-KPO-A	CDLC- 0405-KPO-A	CTLC- 0405-KPO-A	
4	6	CSLC- 0406-KPO-A	CDLC- 0406-KPO-A	CTLC- 0406-KPO-A	
4	8	CSLC- 0408-KPO-A	CDLC- 0408-KPO-A	CTLC- 0408-KPO-A	
4	10	CSLC- 0410-KPO-A	CDLC- 0410-KPO-A	CTLC- 0410-KPO-A	
4	12	CSLC- 0412-KPO-A			
4	13	CSLC- 0413-KPO-A	CDLC- 0413-KPO-A	CTLC- 0413-KPO-A	
4	15	CSLC- 0415-KPO-A	CDLC- 0415-KPO-A		
4	17	CSLC- 0417-KPO-A			
4	20	CSLC- 0420-KPO-A			
4	30			CTLC- 0430-KPO-A	
4.5	8		CDLC- 04.508-KPO-A		
4.5	10		CDLC- 04.510-KPO-A		
5	5		CDLC- 0505-KPO-A	CTLC- 0505-KPO-A	
5	10	CSLC- 0510-KPO-A	CDLC- 0510-KPO-A	CTLC- 0510-KPO-A	
5	13	CSLC- 0513-KPO-A	CDLC- 0513-KPO-A	CTLC- 0513-KPO-A	
5	15	CSLC- 0515-KPO-A	CDLC- 0515-KPO-A	CTLC- 0515-KPO-A	
5	17	CSLC- 0517-KPO-A			
5	20	CSLC- 0520-KPO-A		CTLC- 0520-KPO-A	
5	30			CTLC- 0530-KPO-A	
5.5	10		CDLC- 05.510-KPO-A	CTLC- 05.510-KPO-A	
5.5	15		CDLC- 05.515-KPO-A	CTLC- 05.515-KPO-A	
6	8				CQLC- 0608-KPO-A
6	13				CQLC- 0613-KPO-A
6	15				CQLC- 0615-KPO-A
6	20				CQLC- 0620-KPO-A

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection caps

KIT:CxLC-xxII-KPGY

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen
3	5	CSLC- 0305-KPGY			
3	8	CSLC- 0308-KPGY			
3	13	CSLC- 0313-KPGY			
4	8	CSLC- 0408-KPGY	CDLC- 0408-KPGY	CTLC- 0408-KPGY	
4	13	CSLC- 0413-KPGY	CDLC- 0413-KPGY	CTLC- 0413-KPGY	
5	8		CDLC- 0508-KPGY	CTLC- 0508-KPGY	
5	13		CDLC- 0513-KPGY	CTLC- 0513-KPGY	
5	20			CTLC- 0520-KPGY	
6	8				CQLC- 0608KPGY
6	13				CQLC- 0613KPGY
6	15				CQLC- 0615KPGY
6	20				CQLC- 0620KPGY

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Guiding Y connector
- Fixation Wing

KIT:CxLC-xxII-KPGY-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen	Quadre lumen
3	5	CSLC- 0305-KPGY-A			
3	8	CSLC- 0308-KPGY-A			
3	13	CSLC- 0313-KPGY-A			
4	8	CSLC- 0408-KPGY-A	CDLC- 0408-KPGY-A	CTLC- 0408-KPGY-A	
4	13	CSLC- 0413-KPGY-A	CDLC- 0413-KPGY-A	CTLC- 0413-KPGY-A	
5	8		CDLC- 0508-KPGY-A	CTLC- 0508-KPGY-A	
5	13		CDLC- 0513-KPGY-A	CTLC- 0513-KPGY-A	
5	20			CTLC- 0520-KPGY-A	
6	8				CQLC- 0608KPGY-A
6	13				CQLC- 0613KPGY-A
6	15				CQLC- 0615KPGY-A
6	20				CQLC- 0620KPGY-A

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Guiding Y connector
- Fixation Wing

KIT:CxLC-xxII-KPV

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
2	5	CSLC- 0205-KPV		
2	8	CSLC- 0208-KPV		
3	4		CDLC- 0304-KPV	
3	5	CSLC- 0305-KPV	CDLC- 0305-KPV	
3	8	CSLC- 0308-KPV	CDLC- 0308-KPV	
3	10	CSLC- 0310-KPV	CDLC- 0310-KPV	
3	13	CSLC- 0313-KPV	CDLC- 0313-KPV	
3	15	CSLC- 0315-KPV		
3	17	CSLC- 0317-KPV		
4	5	CSLC- 0405-KPV	CDLC- 0405-KPV	CTLC- 0405-KPV
4	6	CSLC- 0406-KPV	CDLC- 0406-KPV	CTLC- 0406-KPV
4	8	CSLC- 0408-KPV	CDLC- 0408-KPV	CTLC- 0408-KPV
4	10	CSLC- 0410-KPV	CDLC- 0410-KPV	CTLC- 0410-KPV
4	12	CSLC- 0412-KPV		
4	13	CSLC- 0413-KPV	CDLC- 0413-KPV	CTLC- 0413-KPV
4	15	CSLC- 0415-KPV	CDLC- 0415-KPV	
4	17	CSLC- 0417-KPV		
4	20	CSLC- 0420-KPV		
4	30			CTLC- 0430-KPV



KIT:CxLC-xxII-KPV

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
4.5	8		CDLC- 04.508-KPV	
4.5	10		CDLC- 04.510-KPV	
5	5		CDLC- 0505-KPV	CTLC- 0505-KPV
5	6	CSLC- 0506-KPV	CDLC- 0506-KPV	
5	8	CSLC- 0508-KPV	CDLC- 0508-KPV	CTLC- 0508-KPV
5	10	CSLC- 0510-KPV	CDLC- 0510-KPV	CTLC- 0510-KPV
5	13	CSLC- 0513-KPV	CDLC- 0513-KPV	CTLC- 0513-KPV
5	15	CSLC- 0515-KPV	CDLC- 0515-KPV	CTLC- 0515-KPV
5	17	CSLC- 0517-KPV		
5	20	CSLC- 0520-KPV		CTLC- 0520-KPV
5	30			CTLC- 0530-KPV
5.5	10		CDLC- 05.510-KPV	CTLC- 05.510-KPV
5.5	15		CDLC- 05.515-KPV	CTLC- 05.515-KPV

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Dispenser
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection caps

KIT:CxLC-xxII-KPV-A

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
2	5	CSLC- 0205-KPV-A		
2	8	CSLC- 0208-KPV-A		
3	4		CDLC- 0304-KPV-A	
3	5	CSLC- 0305-KPV-A	CDLC- 0305-KPV-A	
3	8	CSLC- 0308-KPV-A	CDLC- 0308-KPV-A	
3	10	CSLC- 0310-KPV-A	CDLC- 0310-KPV-A	
3	13	CSLC- 0313-KPV-A	CDLC- 0313-KPV-A	
3	15	CSLC- 0315-KPV-A		
3	17	CSLC- 0317-KPV-A		
4	5	CSLC- 0405-KPV-A	CDLC- 0405-KPV-A	CTLC- 0405-KPV-A
4	6	CSLC- 0406-KPV-A	CDLC- 0406-KPV-A	CTLC- 0406-KPV-A
4	8	CSLC- 0408-KPV-A	CDLC- 0408-KPV-A	CTLC- 0408-KPV-A
4	10	CSLC- 0410-KPV-A	CDLC- 0410-KPV-A	CTLC- 0410-KPV-A
4	12	CSLC- 0412-KPV-A		
4	13	CSLC- 0413-KPV-A	CDLC- 0413-KPV-A	CTLC- 0413-KPV-A
4	15	CSLC- 0415-KPV-A	CDLC- 0415-KPV-A	
4	17	CSLC- 0417-KPV-A		
4	20	CSLC- 0420-KPV-A		
4	30			CTLC- 0430-KPV-A
4.5	8		CDLC- 04.508-KPV-A	
4.5	10		CDLC- 04.510-KPV-A	
5	5		CDLC- 0505-KPV-A	CTLC- 0505-KPV-A
5	6	CSLC- 0506-KPV-A	CDLC- 0506-KPV-A	
5	8	CSLC- 0508-KPV-A	CDLC- 0508-KPV-A	CTLC- 0508-KPV-A
5	10	CSLC- 0510-KPV-A	CDLC- 0510-KPV-A	CTLC- 0510-KPV-A
5	13	CSLC- 0513-KPV-A	CDLC- 0513-KPV-A	CTLC- 0513-KPV-A
5	15	CSLC- 0515-KPV-A	CDLC- 0515-KPV-A	CTLC- 0515-KPV-A
5	17	CSLC- 0517-KPV-A		
5	20	CSLC- 0520-KPV-A		CTLC- 0520-KPV-A
5	30			CTLC- 0530-KPV-A
5.5	10		CDLC- 05.510-KPV-A	CTLC- 05.510-KPV-A
5.5	15		CDLC- 05.515-KPV-A	CTLC- 05.515-KPV-A

KIT CONTENTS :

- Pediatric Central Venous Catheter
- Introducer needle
- Syringe
- Nitinol guide wire J end with Advancer
- Vessel Dilator
- Scalpel
- Fixation Wing
- Injection caps




4. Umbilical Catheter Kits

KIT:CxLC-xxII-KU


Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
2	40	CSLC-0240-KU		
3	25		CDLC-0325-KU	
3	40	CSLC-0340-KU	CDLC-0340-KU	
4	25		CDLC-0425-KU	CTLC-0425-KU
4	40	CSLC-0440-KU	CDLC-0440-KU	
5	25		CDLC-0525-KU	
5	40	CSLC-0540-KU	CDLC-0540-KU	
8	40	CSLC-0840-KU		

KIT CONTENTS :

- Umbilical Catheter
- Stopcock
- Caps
- Stiffener wire (for sizes 2&3 Fr. Only)



AMECATH



**PERIPHERAL INSERTED
CENTRAL VENOUS CATHETERS
(PICC) KITS**

AMECO MEDICAL INDUSTRIES

Peripheral Inserted Catheter KIT

For direct puncture (KON)

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
3	40	SLPICC -0340-KON		
3	50	SLPICC -0350-KON		
3	55	SLPICC -0355-KON		
3	60	SLPICC -0360-KON		
4	40	SLPICC -0440-KON	DLPICC -0440-KON	
4	50	SLPICC -0450-KON	DLPICC -0450-KON	
4	55	SLPICC -0455-KON	DLPICC -0455-KON	
4	60	SLPICC -0460-KON	DLPICC -0460-KON	
5	40	SLPICC -0540-KON	DLPICC -0540-KON	
5	50	SLPICC -0550-KON	DLPICC -0550-KON	
5	55	SLPICC -0555-KON	DLPICC -0555-KON	
5	60	SLPICC -0560-KON	DLPICC -0560-KON	
6	40		DLPICC -0640-KON	TLPICC -0640-KON
6	50		DLPICC -0650-KON	TLPICC -0650-KON
6	55		DLPICC -0655-KON	TLPICC -0655-KON
6	60		DLPICC -0660-KON	TLPICC -0660-KON

KIT CONTENTS :

- Peripherally inserted Central Catheter
- Touhy-Burst adaptor
- T-connector with sidearm with needless connector
- Nitinol Stiffener with flexible tip
- Over the needle peel away sheath introducer.
- Safety scalpel
- Syringe lock
- Needless connector (valve), 1 for dual lumen and 2 for triple lumen
- Catheter hub securement adhesive. (Unifix lock)
- Transparent catheter adhesives.
- Measuring Tape.


Peripheral inserted Catheter kit

For Seldinger Technique (KOD)

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
3	40	SLPICC -0340-KOD		
3	50	SLPICC -0350-KOD		
3	55	SLPICC -0355-KOD		
3	60	SLPICC -0360-KOD		
4	40	SLPICC -0440-KOD	DLPICC -0440-KOD	
4	50	SLPICC -0450-KOD	DLPICC -0450-KOD	
4	55	SLPICC -0455-KOD	DLPICC -0455-KOD	
4	60	SLPICC -0460-KOD	DLPICC -0460-KOD	
5	40	SLPICC -0540-KOD	DLPICC -0540-KOD	
5	50	SLPICC -0550-KOD	DLPICC -0550-KOD	
5	55	SLPICC -0555-KOD	DLPICC -0555-KOD	
5	60	SLPICC -0560-KOD	DLPICC -0560-KOD	
6	40		DLPICC -0640-KOD	TLPICC -0640-KOD
6	50		DLPICC -0650-KOD	TLPICC -0650-KOD
6	55		DLPICC -0655-KOD	TLPICC -0655-KOD
6	60		DLPICC -0660-KOD	TLPICC -0660-KOD

KIT CONTENTS :

- Peripherally inserted Central Catheter
- Touhy-Burst adaptor
- T-connector with sidearm with needless connector
- Nitinol Stiffener with flexible tip
- Introducer needle with echogenic tip.
- Straight end Nitinol core guidewires
- 5 Cm Over the dilator peel away sheath introducer.
- Safety scalpel
- Syringe lock
- Needleless connector (valve), 1 for dual lumen and 2 for triple lumen
- Catheter hub securement adhesive. (Unifix lock)
- Transparent catheter adhesives.
- Measuring Tape



For Seldinger Technique (KOD)*

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
3	40	SLPICC -0340-KOD*		
3	50	SLPICC -0350-KOD*		
3	55	SLPICC -0355-KOD*		
3	60	SLPICC -0360-KOD*		
4	40	SLPICC -0440-KOD*	DLPICC -0440-KOD*	
4	50	SLPICC -0450-KOD*	DLPICC -0450-KOD*	
4	55	SLPICC -0455-KOD*	DLPICC -0455-KOD*	
4	60	SLPICC -0460-KOD*	DLPICC -0460-KOD*	
5	40	SLPICC -0540-KOD*	DLPICC -0540-KOD*	
5	50	SLPICC -0550-KOD*	DLPICC -0550-KOD*	
5	55	SLPICC -0555-KOD*	DLPICC -0555-KOD*	
5	60	SLPICC -0560-KOD*	DLPICC -0560-KOD*	
6	40		DLPICC -0640-KOD*	TLPICC -0640-KOD*
6	50		DLPICC -0650-KOD*	TLPICC -0650-KOD*
6	55		DLPICC -0655-KOD*	TLPICC -0655-KOD*
6	60		DLPICC -0660-KOD*	TLPICC -0660-KOD*

KIT CONTENTS :

- Peripherally inserted Central Catheter
- Touhy-Burst adaptor
- T-connector with sidearm with needless connector
- 0.018" Nitinol Stiffener with flexible tip
- G21x 7 cm introducer needle with echogenic tip.
- Straight end 0.018" x 50 cm Nitinol core guidewires
- 7 Cm Over the dilator peel away sheath introducer.
- Safety scalpel
- 10 cc Syringe lock
- Needless connector (valve), 1 for dual lumen and 2 for triple lumen
- Catheter hub securement adhesive. (Unifix lock)
- Transparent catheter adhesives.
- Measuring Tape.

2. Pressure Resistant Peripheral Inserted Central Venous Catheter kit

Pressure PICC for Direct Puncture (KON)

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
3	55	PRSLPICC -0355-KON		
4	55	PRSLPICC -0455-KON		
5	55	PRSLPICC -0555-KON	PRDLPICC -0555-KON	
6	55		PRDLPICC -0655-KON	PRTL PICC -0655-KON

KIT CONTENTS :

- Pressure Resistant Peripherally Inserted Central Catheter
- Touhy-Burst adaptor
- T-connector with sidearm with needless connector
- Nitinol Stiffener with flexible tip
- Over the needle peel away sheath introducer.
- Safety scalpel
- Syringe lock
- Needleless connector (valve), 1 for dual lumen and 2 for triple lumen
- Catheter hub securement adhesive. (Unifix lock)
- Transparent catheter adhesives.
- Measuring Tape.

Pressure PICC for Seldinger Technique (KOD)

Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
3	55	PRSLPICC -0355-KOD		
4	55	PRSLPICC -0455-KOD		
5	55	PRSLPICC -0555-KOD	PRDLPICC -0555-KOD	
6	55		PRDLPICC -0655-KOD	PRTL PICC -0655-KOD

KIT CONTENTS :


- Pressure resistant Peripherally inserted Central Catheter
- Touhy-Burst adaptor
- T-connector with sidearm with needless connector
- Nitinol Stiffener with flexible tip
- Introducer needle with echogenic tip.
- Straight end Nitinol core guidewires
- 5 Cm Over the dilator peel away sheath introducer.
- Safety scalpel
- Syringe lock
- Needleless connector (valve), 1 for dual lumen and 2 for triple lumen
- Catheter hub securement adhesive. (Unifix lock)
- Transparent catheter adhesives.
- Measuring Tape.

Pressure PICC for Seldinger Technique (KOD)*


Fr. Size	Length	Single lumen	Dual lumen	Triple lumen
3	55	PRSLPICC -0355-KOD*		
4	55	PRSLPICC -0455-KOD*		
5	55	PRSLPICC -0555-KOD*	PRDLPICC -0555-KOD*	
6	55		PRDLPICC -0655-KOD*	PRTL PICC -0655-KOD*

KIT CONTENTS :

- Pressure resistant Peripherally inserted Central Catheter
- Touhy-Burst adaptor
- T-connector with sidearm with needless connector
- Nitinol Stiffener with flexible tip
- Introducer needle with echogenic tip.
- Straight end Nitinol core guidewires
- 7 Cm Over the dilator peel away sheath introducer.
- Safety scalpel
- Syringe lock
- Needleless connector (valve), 1 for dual lumen and 2 for triple lumen
- Catheter hub securement adhesive. (Unifix lock)
- Transparent catheter adhesives.
- Measuring Tape.



AMECATH



MIDLINE VENOUS CATHETER KIT

AMECO MEDICAL INDUSTRIES

1. Midline Catheter Kit for Seldinger Technique (using IV Cath Cannula) :

Fr. Size	Length (cm)	Single Lumen	Dual Lumen	Pressure Resistant Single Lumen	Pressure Resistant Dual Lumen
3	10	MLSLC-0310-K		PRMLSLC-0310-K	
3	12	MLSLC-0312-K		PRMLSLC-0312-K	
3	20	MLSLC-0320-K		PRMLSLC-0320-K	
3	25	MLSLC-0325-K		PRMLSLC-0325-K	
4	10	MLSLC-0410-K		PRMLSLC-0410-K	
4	12	MLSLC-0412-K		PRMLSLC-0412-K	
4	20	MLSLC-0420-K	MLDLC-0420-K	PRMLSLC-0420-K	PRMLDLC-0420-K
4	25	MLSLC-0425-K	MLDLC-0425-K	PRMLSLC-0425-K	PRMLDLC-0425-K
5	12	MLSLC-0512-K		PRMLSLC -0512-K	
5	20	MLSLC-0520-K	MLDLC-0520-K	PRMLSLC-0520-K	PRMLDLC-0520-K
5	25	MLSLC-0525-K	MLDLC-0525-K	PRMLSLC-0525-K	PRMLDLC-0525-K

KIT CONTENTS :

- Midline Venous Catheter
- Over the Needle Catheter
- Straight end Nitinol Guide Wire
- Needless Connector (one for single lumen and two for dual lumen)



2. Midline Catheter Kit for Seldinger Technique (using over needle peel away introducer) :

Fr. Size	Length (cm)	Single Lumen	Dual Lumen	Pressure Resistant Single Lumen	Pressure Resistant Dual Lumen
3	10	MLSLC-0310-KON		PRMLSLC-0310-KON	
3	12	MLSLC-0312-KON		PRMLSLC-0312-KON	
3	20	MLSLC-0320-KON		PRMLSLC-0320-KON	
3	25	MLSLC-0325-KON		PRMLSLC-0325-KON	
4	10	MLSLC-0410-KON		PRMLSLC-0410-KON	
4	12	MLSLC-0412-KON		PRMLSLC-0412-KON	
4	20	MLSLC-0420-KON	MLDLC-0420-KON	PRMLSLC-0420-KON	PRMLDLC-0420-KON
4	25	MLSLC-0425-KON	MLDLC-0425-KON	PRMLSLC-0425-KON	PRMLDLC-0425-KON
5	12	MLSLC-0512-KON		PRMLSLC-0512-KON	
5	20	MLSLC-0520-KON	MLDLC-0520-KON	PRMLSLC-0520-KON	PRMLDLC-0520-KON
5	25	MLSLC-0525-KON	MLDLC-0525-KON	PRMLSLC-0525-KON	PRMLDLC-0525-KON

KIT CONTENTS :

- Midline Venous Catheter
- 0.018” Nitinol Stiffener with Flexible end
- Over the Needle Peel Away Sheath Introducer.
- Touhy Burst Adapter
- T-Connector with Side Arm with Needless Connector
- 10 cc Syringe
- Safety Scalpel
- Needless Connector (with dual lumen catheters only)
- Catheter Hub Securement Adhesive. (Unifix Lock)
- Transparent Catheter Adhesive.

3. Midline Catheter Kit for Seldinger Technique (using over dilator peel away introducer) :

Fr. Size	Length	Single Lumen	Dual Lumen	Pressure Resistant Single Lumen	Pressure Resistant Dual Lumen
3	10	MLSLC-0310-KOD		PRMLSLC-0310-KOD	
3	12	MLSLC-0312-KOD		PRMLSLC-0312-KOD	
3	20	MLSLC-0320-KOD		PRMLSLC-0320-KOD	
3	25	MLSLC-0325-KOD		PRMLSLC-0325-KOD	
4	10	MLSLC-0410-KOD		PRMLSLC-0410-KOD	
4	12	MLSLC-0412-KOD		PRMLSLC-0412-KOD	
4	20	MLSLC-0420-KOD	MLDLC-0420-KOD	PRMLSLC-0420-KOD	PRMLDLC-0420-KOD
4	25	MLSLC-0425-KOD	MLDLC-0425-KOD	PRMLSLC-0425-KOD	PRMLDLC-0425-KOD
5	12	MLSLC-0512-KOD		PRMLSLC-0512-KOD	
5	20	MLSLC-0520-KOD	MLDLC-0520-KOD	PRMLSLC-0520-KOD	PRMLDLC-0520-KOD
5	25	MLSLC-0525-KOD	MLDLC-0525-KOD	PRMLSLC-0525-KOD	PRMLDLC-0525-KOD

KIT CONTENTS :

- Midline Venous Catheter
- 0.018" Nitinol Stiffener with Flexible Tip
- G21 x 7 cm Introducer Needle with Echogenic Tip
- Touhy Burst Adapter
- T-Connector with Side Arm with Needless Connector
- 10 cc Syringe
- Straight end 0.018"x 50 cm Nitinol Guidewire.
- Over the Dilator Peel Away Sheath Introducer.
- Safety Scalpel
- Needless Connector (with dual lumen catheters only)
- Catheter Hub Securement Adhesive. (Unifix Lock)
- Transparent Catheter Adhesive.

AMECATH



ARTERIAL CATHETER KIT

AMECO MEDICAL INDUSTRIES

1. Arterial catheters. For Seldinger Technique kit (KS) :

Size (Fr)	Size (Ga)	Length (cm)	CODE
3	20	5	AL-0305-KS
3	20	6	AL-0306-KS
3	20	8	AL-0308-KS
3	20	12	AL-0312-KS
4	18	8	AL-0408-KS
4	18	10	AL-0410-KS
4	18	16	AL-0416-KS
5	16	8	AL-0508-KS
5	16	12	AL-0512-KS
5	16	16	AL-0516-KS

KIT CONTENTS :

- Arterial Catheter
- Introducer Needle
- Straight end Nitinol Guide Wire
- Cap

2. Arterial catheters. For Seldinger Technique kit with Disposable Cartridge (KSD) :

Size (Fr)	Size (Ga)	Length (cm)	CODE
3	20	5	AL-0305-KSD
3	20	6	AL-0306-KSD
3	20	8	AL-0308-KSD
3	20	12	AL-0312-KSD
4	18	8	AL-0408-KSD
4	18	10	AL-0410-KSD
4	18	16	AL-0416-KSD
5	16	8	AL-0508-KSD
5	16	12	AL-0512-KSD
5	16	16	AL-0516-KSD

KIT CONTENTS :

- Arterial Catheter
- Introducer Needle
- Straight end Nitinol Guide Wire
- Disposable Cartridge
- Cap

- Established at 2001 for design, development and production of sterile dialysis catheter kits and accessories, Central venous kits and accessories, Urological catheters kits, stents, needles and accessories.
- 250 people are currently working at Ameco among them 5 doctors, 8 engineers, high caliber administrative, logistic and production .
- All manufacturing processes are under one roof. Extrusion, molding, over molding, tipping, punching, assembly, coating packaging and sterilization.
- Many Success stories for contact manufacturing and private labeling with major international companies.
- In 2002, Amecath received ISO13485 certificate in scope of design, development and production of catheter and kits, entailing that Amecath is applying full product quality assurance and Quality management system from design, production to service provision.
- In 2002, Amecath received EC mark (EC certificate) for all products, entailing that Amecath products can be freely placed in the market of 25 European countries and other countries with recognition of the European directives
- In 2005, Amecath received EC-Design examination for class III products (Central Venous Haemodialysis catheters for long term use), entailing that the device has passed the examination for safety and performance.
- In 2011, Amecath received Health Canada approval for central venous catheters, short term haemodialysis catheter and sheath introducer and urology double loop stents.
- In 2011, Amecath received Russian registration for central venous catheter and hemodialysis catheter, other products are under the registration process
- In 2012, Amecath received FDA approval for central venous catheter and short term hemodialysis catheter other products are currently under fda review.
- In 2013, Amecath received FDA approval for Dual lumen implantable long term haemodialysis catheters, sheath introducers and other products are currently under FDA review.
- In 2013, Amecath received EC-Design examination for class III products (Central Venous and peripherally inserted central catheters), entailing that the device has passed the examination for safety and performance
- In 2013, Amecath received health Canada Approval for long term implantable hemodialysis catheters

Manufactured By:**AMECO MEDICAL INDUSTRIES**

Industrial Zone B4 Plot 119 east,
10th of Ramadan city - Egypt

Tel: +20.554.501321/ 2

Fax :+20.554.501224

Mail: support@amecathgroup.com

Website: www.amecathgroup.com

Authoized European Representative**Obelis s.a**

Bd. Général Wahis 53
1030 Brussels, BELGIUM

Tel: +32.2. 732.59.54

Fax: +32.2.732.60.03

Email: mail@obelis.net

Website: www.obelis.net