

# Pharmapad & Pharmapad Carbon



**Pharmapad** is a wound dressing, having a low-adherent finely perforated net on both sides of an absorbent wound pad.

**Pharmapad Carbon** is a low adherent multi-layered absorbent wound pad with additional non-woven layer containing activated carbon for managing wounds with offensive odour.

## Features:

- Absorbent
- Low-adherent wound contact layer.
- Air permeable and conformable.
- The perforated wound contact layer allows rapid passage of wounds exudates to the absorbent core.
- Product can be used on both sides, so mistakes of application is impossible.
- The additional activated carbon layer effectively adsorb wound's offensive odors.
- Can be cut to the required size.

## Intended purpose:

**Pharmapad family** is indicated for the protection and absorption of wound exudates from:

- Lightly exuding wounds.
- Clean post-operative incisions.
- For first aid use, e.g. cuts and abrasions.
- Second degree burns.
- Variants containing carbon are indicated for malodorous acute or chronic wounds.

# Pharmapad & Pharmapad Carbon

Pharmapad	Pharmapad Carbon	size	pieces / box	Boxes/Case	Case Dimensions
PAD5050	PADC5050	5cmx5cm	100	16	67.5x30x33 cm
PAD5075	PADC5075	5cmx7.5cm	100	16	67.5x30x33 cm
PAD75100	PADC75100	7.5cmx10cm	50	16	67.5x30x33 cm
PAD75150	PADC75150	7.5cmx15cm	50	16	67.5x30x33 cm
PAD75200	PADC75200	7.5cmx20cm	40	16	67.5x30x33 cm
PAD100100	PADC100100	10cmx10cm	50	16	67.5x30x33 cm
PAD100150	PADC100150	10cmx15cm	50	16	67.5x30x33 cm
PAD100200	PADC100200	10cmx20cm	40	16	67.5x30x33 cm
PAD50100	PADC50100	5cmx10cm	100	16	67.5x30x33 cm
PAD100250	PADC100250	10cmx25cm	40	16	67.5x30x33 cm
PAD100300	PADC100300	10cmx30cm	40	16	67.5x30x33 cm
PAD150250	PADC150250	15cmx25cm	10	24	53x27x27 cm
PAD20300	PADC20300	2cmx30cm	5	24	47x47x20 cm
PAD20400	PADC20400	2cmx40cm	5	24	47x47x20 cm
PAD200200	PADC200200	20cmx20cm	10	24	53x27x27 cm
PAD200300	PADC200300	20cmx30cm	10	24	53x27x27 cm
PAD25450	PADC25450	2.5cmx45cm	5	24	47x47x20 cm
PAD50300	PADC50300	5cmx30cm	50	16	67.5x30x33 cm
PAD70120	PADC70120	7cmx12cm	50	16	67.5x30x33 cm
PAD75120	PADC75120	7.5cmx12cm	50	16	67.5x30x33 cm
PAD7575	PADC7575	7.5cmx7.5cm	50	16	67.5x30x33 cm
PAD150150	PADC150150	15cmx15cm	10	24	53x27x27 cm
PAD150200	PADC150200	15cmx20cm	10	24	53x27x27 cm
PAD200400	PADC200400	20x40cm	5	24	67.5 x 30 x 33 cm
PAD5070	PADC5070	5x7 cm	100	16	67.5 x 30 x 33 cm
PAD6565T		6.5x6.5 Tracheostomy	50	16	67.5 x 30 x 33 cm
PAD6080T		6x8 cm Tracheostomy	50	16	67.5 x 30 x 33 cm
PAD8090T		8x9 cm Tracheostomy	50	16	67.5 x 30 x 33 cm
PAD5882		5.8 x 8.2 cm	100	16	

- Amria - Alexandria - Egypt  
 Tel. : 002 03 4500264 Fax : 002 03 4500263  
 - Kafr-El-Zayat P.O.Box 38 Egypt  
 Tel. : 002 040 2545643 Fax : 002 040 2545645  
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**  
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3 )						
Customer : Pharmaplast		Product name : pharmaapd & phramapad carbon						
Issue date : 13/06/2015		Prepared by : Nader Nabil Zaki						
Colors : <input type="checkbox"/>		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	02/05/2019	new text by jostina	(Essential requirement to be updated with clinical evaluation update					
3	29-08-2021	Addition of the new codes	To be aligned with the new coding system					
4	18-05-2022	table of sizes updates	File maintenance					