

PrecisionCut Ophthalmic MicroSurgical Blades (Small Incision Surgical Blades)

- ▶ Made from medical grade stainless steel (SS420)
- ▶ Precision craftsmanship with accuracy of angle, thickness, bends and sizes
- ▶ Ultra precision sharpness
- ▶ Supplied in a special storage case so as to prevent accidental damage to tips
- ▶ Distinct colored handles for easy identification
- ▶ Meets distinct requirements for every ophthalmic surgery with wide spectrum of selection to choose from
- ▶ The Action PrecisionCut Small Incision Surgical Blade range offers an optimal combination of Quality & Economy
We offer all variants (sizes) of blades in Keratome, Crescent, Side Port (Lance Tip), Clear Cornea, MVR Blades
- ▶ Single Bevel, Double Bevel, Clear Cornea, Depth Control



- ▶ Customize your own colour handles

