

# Protectfilm® Family



## Description:

**Protectfilm® family** is a transparent protective wound dressing that consists of thin bacteria-proof, virus-proof and waterproof polyurethane film coated with hypoallergenic adhesive.

Available as **Protectfilm® Frame Style** delivery system for easier handling.

**Protectfilm® Adhesive border** is a highly breathable transparent protective polyurethane film with adhesive border and adhesive free central area which prevents undesired adhesion with the wound bed and allows removal of the dressing without pain or trauma.

## Features:

- Waterproof: **Protectfilm®** allows the patient to bathe without dressing change.
- Comfortable: as it conforms to body contours and is easy to apply.
- Transparent: allows wound monitoring without disturbing the injury or removing the dressing.
- Bacteria-proof and virus- proof: protects the wound from external contaminants.
- Breathable: **Protectfilm® Adhesive border** with its adhesive free central area, presents a much higher breathability than most ordinary wound dressings.

## Intended Purpose:

- As a primary dressing, it's intended to cover and protect dry or slightly exuding wounds such as: abrasions, lacerations and areas of friction.
- As a secondary dressing, it's used as a fixation device, to secure catheter sites and to cover other primary dressings like alginates, hydrogels, hydrofibers and foams.

# Protectfilm® Family

Protectfilm®	Size	Pieces/Box	Boxes/Case	Case Dimensions
PROT6070	6cmx7cm	100	16	67.5X30X33 cm
PROT5075	5cmx7.5cm	100	16	67.5X30X33 cm
PROT75100	7.5cmx10cm	100	16	67.5X30X33 cm
PROT75120	7.5cmx12cm	60	16	67.5X30X33 cm
PROT5082	5cmx8.2cm	100	16	67.5X30X33 cm
PROT6585	6.5cmx8.5cm	100	16	67.5X30X33 cm
PROT100120	10cmx12cm	60	16	67.5X30X33 cm
PROT100140	10cmx14cm	60	16	67.5X30X33 cm
PROT100160	10cmx16cm	60	16	67.5X30X33 cm
PROT4545	4.5cmx4.5cm	100	16	67.5X30X33 cm
PROT4867	4.8cmx6.7cm	100	16	67.5X30X33 cm
PROT100100	10cmx10cm	100	16	67.5X30X33 cm
PROT100150	10cmx15cm	60	16	67.5X30X33 cm
PROT100200	10cmx20cm	60	16	67.5X30X33 cm
PROT100250	10cmx25cm	60	16	67.5X30X33 cm
PROT100300	10cmx30cm	60	8	44X30X30 cm
PROT150200	15cmx20cm	50	16	67.5X30X33 cm
PROT200300	20cmx30cm	25	8	44X30X30 cm
PROT150150	15cmx15cm	25	8	44X30X30 cm
PROT200200	20cmx20cm	25	8	44X30X30 cm
PROT165180	16.5cmx18cm	25	8	44X30X30 cm
PROT220220	22cmx22cm	25	8	44X30X30 cm
PROT100180	10cmx18cm	50	16	67.5X30X33 cm
PROT150250	15cmx25cm	25	8	44X30X30 cm
PROT200250	20cmx25cm	25	16	67.5X30X33 cm
PROT4870	4.8cmx7cm	100	16	67.5X30X33 cm
PROT5070	5cmx7cm	100	16	67.5X30X33 cm
PROT5072	5cmx7.2cm	100	16	67.5X30X33 cm
PROT5080	5cmx8cm	100	16	67.5X30X33 cm
PROT6080	6cmx8cm	100	16	67.5X30X33 cm

## Protectfilm® Family

Protectfilm® Frame Style	Size	Pieces/Box	Boxes/Case	Case Dimensions
PRFS6070	6cmx7cm	100	16	58X32X13 cm
PRFS5075	5cmx7.5cm	100	16	58X32X13 cm
PRFS75100	7.5cmx10cm	100	16	44X30X30 cm
PRFS75120	7.5cmx12cm	60	16	44X30X30 cm
PRFS5082	5cmx8.2cm	100	16	58X32X13 cm
PRFS6585	6.5cmx8.5cm	100	16	58X32X13 cm
PRFS100120	10cmx12cm	60	16	44X42X17 cm
PRFS100140	10cmx14cm	60	16	44X30X30 cm
PRFS100160	10cmx16cm	60	16	44X30X30 cm
PRFS4545	4.5cmx4.5cm	100	16	58X32X13 cm
PRFS4867	4.8cmx6.7cm	100	16	58X32X13 cm
PRFS100100	10cmx10cm	100	16	44X30X30 cm
PRFS100150	10cmx15cm	60	16	44X30X30 cm
PRFS100200	10cmx20cm	60	16	44X30X30 cm
PRFS100250	10cmx25cm	60	16	44X30X30 cm
PRFS100300	10cmx30cm	60	8	
PRFS150200	15cmx20cm	50	16	
PRFS200300	20cmx30cm	25	8	
PRFS150150	15cmx15cm	25	8	
PRFS200200	20cmx20cm	25	8	
PRFS165180	16.5cmx18cm	25	8	
PRFS220220	22cmx22cm	25	8	
PRFS100180	10cmx18cm	50	16	44X30X30 cm
PRFS150250	15cmx25cm	25	8	
PRFS200250	20cmx25cm	25	16	
PRFS4870	4.8cmx7cm	100	16	58X32X13 cm
PRFS5070	5cmx7cm	100	16	58X32X13 cm
PRFS5072	5cmx7.2cm	100	16	58X32X13 cm
PRFS5080	5cmx8cm	100	16	58X32X13 cm
PRFS6080	6cmx8cm	100	16	58X32X13 cm

# Protectfilm® Family

Protectfilm® Adhesive Border	Size	Pieces/Box	Boxes/Case	Case Dimensions
PFZ6070ADH	6cmx7cm	100	16	67.5x30x32 cm
PFZ5075ADH	5cmx7.5cm	100	16	67.5x30x32 cm
PFZ75100ADH	7.5cmx10cm	100	16	67.5x30x32 cm
PFZ75120ADH	7.5cmx12cm	60	16	67.5x30x32 cm
PFZ5082ADH	5cmx8.2cm	100	16	67.5x30x32 cm
PFZ6585ADH	6.5cmx8.5cm	100	16	67.5x30x32 cm
PFZ100120ADH	10cmx12cm	60	16	67.5x30x32 cm
PFZ100140ADH	10cmx14cm	60	16	67.5x30x32 cm
PFZ100160ADH	10cmx16cm	60	16	67.5x30x32 cm
PFZ4545ADH	4.5cmx4.5cm	100	16	67.5x30x32 cm
PFZ4867ADH	4.8cmx6.7cm	100	16	67.5x30x32 cm
PFZ100100ADH	10cmx10cm	100	16	67.5x30x32 cm
PFZ100150ADH	10cmx15cm	60	16	67.5x30x32 cm
PFZ100200ADH	10cmx20cm	60	16	67.5x30x32 cm
PFZ100250ADH	10cmx25cm	60	16	67.5x30x32 cm
PFZ100300ADH	10cmx30cm	60	8	67.5x30x32 cm
PFZ150200ADH	15cmx20cm	50	16	67.5x30x32 cm
PFZ200300ADH	20cmx30cm	25	8	67.5x30x32 cm
PFZ150150ADH	15cmx15cm	25	8	67.5x30x32 cm
PFZ200200ADH	20cmx20cm	25	8	67.5x30x32 cm
PFZ165180ADH	16.5cmx18cm	25	8	67.5x30x32 cm
PFZ220220ADH	22cmx22cm	25	8	67.5x30x32 cm
PFZ100180ADH	10cmx18cm	50	16	67.5x30x32 cm
PFZ150250ADH	15cmx25cm	25	8	67.5x30x32 cm
PFZ200250ADH	20cmx25cm	25	16	67.5x30x32 cm
PFZ4870ADH	4.8cmx7cm	100	16	67.5x30x32 cm
PFZ5070ADH	5cmx7cm	100	16	67.5x30x32 cm
PFZ5072ADH	5cmx7.2cm	100	16	67.5x30x32 cm
PFZ5080ADH	5cmx8cm	100	16	67.5x30x32 cm
PFZ6080ADH	6cmx8cm	100	16	67.5x30x32 cm

- Amria - Alexandria - Egypt  
 Tel. : 002 03 4500264 Fax : 002 03 4500263  
 - Kafr-El-Zayat P.O.Box 38 Egypt  
 Tel. : 002 040 2545643 Fax : 002 040 2545645  
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**  
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : protectfilm family						
Issue date : 31/08/2016		Prepared by : Marina Medhat						
Colors : <input type="checkbox"/>		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	11/12/2018	We did amendments in packaging capacity	Essential requirement					
3	29/01/2020	amend text and table of size	file maintenance					
4	27/02/2020	amend table of size	add carton size					
5	06-05-2020	Adding more sizes to table of	file maintenance					
6	18-05-2020	adding new packaging configuration according to CRS	Customer request					
7	30-07-2022	- Updating codes - Adding case dimensions	To be aligned with the new coding system					