

AMECATH



DIALYSIS CATHETERS

AMECO MEDICAL INDUSTRIES



DIALYSIS

Catheters and Accessories

1. Short Term Haemodialysis Catheters	2
1.1. Single Lumen Catheter - Desilet Kit	3
1.2. Double Lumen catheters	4
1.2.1. Tapered Tip Dual Lumen Catheter Kits	4
• STD kit	4
• Guiding Syringe kit	5
1.2.2. UniQath	6
• STD kit	6
• Guiding Syringe kit	6
1.2.3. Soft Step Tip Soft Dual Lumen Catheter Kits	7
• STD kit	7
• Guiding Syringe kit	8
1.2.4. Smart Dual Lumen Haemodialysis Catheter Kits	9
• STD kit	9
• Guiding Syringe kit	9
1.3. Triple Lumen catheters	10
1.3.1. Triple Lumen catheter kits	11
• STD kit	11
• Guiding Syringe kit	11
1.3.2. Smart Triple Lumen catheter kits	12
• STD kit	12
• Guiding Syringe kit	12
2. Long Term haemodialysis Catheters	13
2.1. Dual Lumen catheters	14
2.2. Y-Split catheters	15
2.3. Twin Cath catheters	16
2.3.1 P2TC with Integrated Female Luer Valve Kit	16
2.3.2 Twin lumen Haemodialysis Catheter Kit (Pediatric)	16
2.3.3 Twin lumen Haemodialysis Catheter Kit (Adult)	17
2.4 X- Split catheters	18
3. Peritoneal Catheters	19
3.1. Permanent Straight Catheter	20
3.2. Permanent Curled Catheter	21

1. Short Term Haemodialysis Catheter

An Optimum Compromise between competing Factors. Where a balance between:

- Flexibility & firmness.
- Softness and anti kinking.
- Material and durability.
- Size and flow.

Features and Properties

CATHETER FEATURES:

Soft tip.

Length marked.

Small rotating fixation wing.

Extension lines stand repeated clamping and declamping. Multiple threads female luer end.

Injection caps.

Variety of sizes, lengths and shapes.


KIT FEATURES:

Heavy duty length marked guide wire with advancer.

Long tapered tip dilators with stopper to avoid over insertion. Scalpel # 11 for small snipping.

Guiding Y Connector "GY" (option).

Nitinol Guide wire "N".



1.1. Single Lumen Catheter Desilet Kit

Ordering Information:

Fr Size	Length	Desilet
10	13 Cm	SLC-1013-DK
10	15 Cm	SLC-1015-DK
10	18 Cm	SLC-1018-DK
10	20 Cm	SLC-1020-DK

Kit Contents

- (1) Polyurethane Single lumen Haemodialysis catheter.
- (1) J End Nitinol guide wire in scaled dispenser.
- (1) Vessel dilator
- (1) Syringe
- (1) Introducer needle, echogenic
- (1) Scalpel
- (1) Injection cap
- (1) Towel



1.2. Double Lumen catheters

1.2.1 Tapered Tip Dual Lumen Catheter Kits

Ordering Information:
STD Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
6	08Cm	DLC-0608-KP		
7	10Cm	DLC-0710-KP	DLC-0710-KJP	
7	13Cm	DLC-0713-KP	DLC-0713-KJP	
8	10Cm	DLC-0810-KP	DLC-0810-KJP	
8	12Cm	DLC-0812-KP	DLC-0812-KJP	
8	13Cm	DLC-0813-KP	DLC-0813-KJP	
9	10Cm	DLC-0910-KP	DLC-0910-KJP	
9	12Cm	DLC-0912-KP	DLC-0912-KJP	
9	13Cm	DLC-0913-KP	DLC-0913-KJP	
9	13Cm	DLC-0913-K	DLC-0913-KJ	DLC-0913-KPC
9	15Cm	DLC-0915-K	DLC-0915-KJ	DLC-0915-KPC
10	10Cm	DLC-1010-K	DLC-1010-KJ	
10	13Cm	DLC-1013-K	DLC-1013-KJ	DLC-1013-KPC
10	15Cm	DLC-1015-K	DLC-1015-KJ	DLC-1015-KPC
12	12Cm	DLC-1212-K	DLC-1212-KJ	DLC-1212-KPC
12	13Cm	DLC-1213-K	DLC-1213-KJ	DLC-1213-KPC
12	15Cm	DLC-1215-K	DLC-1215-KJ	DLC-1215-KPC
12	16Cm	DLC-1216-K	DLC-1216-KJ	DLC-1216-KPC
12	17Cm	DLC-1217-K	DLC-1217-KJ	DLC-1217-KPC
12	20Cm	DLC-1220-K	DLC-1220-KJ	DLC-1220-KPC
12	24Cm	DLC-1224-K	DLC-1224-KJ	DLC-1224-KPC
14	15Cm	DLC-1415-K	DLC-1415-KJ	DLC-1415-KPC
14	17Cm	DLC-1417-K	DLC-1417-KJ	DLC-1417-KPC
14	20Cm	DLC-1420-K	DLC-1420-KJ	DLC-1420-KPC
16	15Cm	DLC-1615-K	DLC-1615-KJ	
16	17Cm	DLC-1617-K	DLC-1617-KJ	
16	20Cm	DLC-1620-K	DLC-1620-KJ	

Kit Contents

- (1) Tapered Tip Dual Lumen Short Term Haemodialysis Catheter.
- (1) J End Nitinol Guide Wire with Scaled Dispenser
- (2) Vessel Dilators
- (1) Syringe
- (1) Guiding Y Connector
- (1) Introducer Needle, Echogenic
- (1) Scalpel
- (2) Injection Cap
- (1) Movable Fixation Wing (with Pre-curved catheter only)



Maximum Flow - Maximum comfort

- Soft polyurethane for patient comfort
- Stepped tip for maximum flow.

Ordering Information: GS Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
9	13Cm	DLC-0913-K-GS	DLC-0913-KJ-GS	DLC-0913-KPC-GS
9	15Cm	DLC-0915-K-GS	DLC-0915-KJ-GS	DLC-0915-KPC-GS
10	10Cm	DLC-1010-K-GS	DLC-1010-KJ-GS	
10	13Cm	DLC-1013-K-GS	DLC-1013-KJ-GS	DLC-1013-KPC-GS
10	15Cm	DLC-1015-K-GS	DLC-1015-KJ-GS	DLC-1015-KPC-GS
12	12Cm	DLC-1212-K-GS	DLC-1212-KJ-GS	DLC-1212-KPC-GS
12	13Cm	DLC-1213-K-GS	DLC-1213-KJ-GS	DLC-1213-KPC-GS
12	15Cm	DLC-1215-K-GS	DLC-1215-KJ-GS	DLC-1215-KPC-GS
12	16Cm	DLC-1216-K-GS	DLC-1216-KJ-GS	DLC-1216-KPC-GS
12	17Cm	DLC-1217-K-GS	DLC-1217-KJ-GS	DLC-1217-KPC-GS
12	20Cm	DLC-1220-K-GS	DLC-1220-KJ-GS	DLC-1220-KPC-GS
12	24Cm	DLC-1224-K-GS	DLC-1224-KJ-GS	DLC-1224-KPC-GS
14	15Cm	DLC-1415-K-GS	DLC-1415-KJ-GS	DLC-1415-KPC-GS
14	17Cm	DLC-1417-K-GS	DLC-1417-KJ-GS	DLC-1417-KPC-GS
14	20Cm	DLC-1420-K-GS	DLC-1420-KJ-GS	DLC-1420-KPC-GS
16	15Cm	DLC-1615-K-GS	DLC-1615-KJ-GS	
16	17Cm	DLC-1617-K-GS	DLC-1617-KJ-GS	
16	20Cm	DLC-1620-K-GS	DLC-1620-KJ-GS	

Kit Contents

- (1) Tapered Tip Dual Lumen Short Term Haemodialysis Catheter.
- (1) J End Nitinol Guide Wire with Scaled Dispenser
- (2) Vessel Dilators
- (1) Guiding Syringe
- (1) Introducer Needle, Echogenic
- (1) Scalpel
- (2) Injection Cap
- (1) Movable Fixation Wing (with Pre-curved catheter only)



1.2.2. UniQath

New Soft short to midterm haemodialysis catheter. One catheter suitable for all approaches whether Jugular, Subclavian or femoral. The catheter offers high patient comfort, entire haemostatic insertion, months of use, low incidence of infection and reliable flow rate. Being with a fixed length and size, it is one catheter that fit with all patients.

Features and Properties

CATHETER FEATURES:

Soft material: comfortable for the patient & easily inserted. Smart groove: protect side holes from sticking to vessel wall. Tubular design: No turbulence, higher flow.

No connection: prevent turbulence, more comfort.

No clamp: long term saved integrity of extensions.

Haemostatic valves: prevent blood leakage all over the procedure.

Ordering Information: STD Kit

Fr. Size	Length	Subclavian Kit Code
12	15 Cm	UN-DLC-1215-K
12	17 Cm	UN-DLC-1217-K
12	20 Cm	UN-DLC-1220-K
12	25 Cm	UN-DLC-1225-K
14	15 Cm	UN-DLC-1415-K
14	17 Cm	UN-DLC-1417-K
14	20 Cm	UN-DLC-1420-K
14	25 Cm	UN-DLC-1425-K

Ordering Information: GS Kit

Fr. Size	Length	Subclavian Kit Code
12	15 Cm	UN-DLC-1215-K-GS
12	17 Cm	UN-DLC-1217-K-GS
12	20 Cm	UN-DLC-1220-K-GS
12	25 Cm	UN-DLC-1225-K-GS
14	15 Cm	UN-DLC-1415-K-GS
14	17 Cm	UN-DLC-1417-K-GS
14	20 Cm	UN-DLC-1420-K-GS
14	25 Cm	UN-DLC-1425-K-GS



Kit Contents

- (1) Dual Lumen Haemodialysis Catheter
- Stiffener with Haemostatic Valve
- (1) J End Nitinol Guide Wire with Scaled Dispenser
- (2) Vessel Dilators
- (1) Syringe
- (1) Guiding Y Connector
- (1) Introducer Needle, Echogenic
- (1) Scalpel
- (2) Injection Cap
- (1) Unifix Catheter Tube Fixation Adhesive
- (2) Transparent Catheter Fixation Adhesive
- (1) Movable Fixation Wing

Kit Contents

- (1) Dual Lumen Haemodialysis Catheter
- Stiffener with Haemostatic Valve
- (1) J End Nitinol Guide Wire with Scaled Dispenser
- (2) Vessel Dilators
- (1) Guiding Syringe
- (1) Introducer Needle, Echogenic
- (1) Scalpel
- (2) Injection Cap
- (1) Unifix Catheter Tube Fixation Adhesive
- (2) Transparent Catheter Fixation Adhesive
- (1) Movable Fixation Wing


1.2.3 Soft Step Tip Soft Dual Lumen Catheter Kits

Ordering Information:
STD Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
12	12Cm	SDLC-1212-K	SDLC-1212-KJ	
12	15Cm	SDLC-1215-K	SDLC-1215-KJ	SDLC-1215-KPC
12	17Cm	SDLC-1217-K	SDLC-1217-KJ	SDLC-1217-KPC
12	20Cm	SDLC-1220-K	SDLC-1220-KJ	SDLC-1220-KPC
12	25Cm	SDLC-1225-K	SDLC-1225-KJ	SDLC-1225-KPC
14	15Cm	SDLC-1415-K	SDLC-1415-KJ	SDLC-1415-KPC
14	17Cm	SDLC-1417-K	SDLC-1417-KJ	SDLC-1417-KPC
14	20Cm	SDLC-1420-K	SDLC-1420-KJ	SDLC-1420-KPC
14	25Cm	SDLC-1425-K	SDLC-1425-KJ	SDLC-1425-KPC
14	30Cm	SDLC-1430-K	SDLC-1430-KJ	SDLC-1430-KPC
14	35Cm	SDLC-1435-K	SDLC-1435-KJ	SDLC-1435-KPC

Kit Contents

- (1) Step Tip Soft Dual Lumen Haemodialysis Subclavian Catheter
- Stiffener with Haemostatic Valve
- (1) J End Nitinol Guide Wire with Scaled Dispenser
- (2) Vessel Dilators
- (1) Syringe
- (1) Guiding Y Connector
- (1) Introducer Needle, Echogenic
- (1) Scalpel
- (2) Injection Cap
- (1) Unifix Catheter Tube Fixation Adhesive
- (2) Transparent Catheter Fixation Adhesive
- (1) Movable Fixation Wings (with Pre-curved catheter only)



Ordering Information: GS Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
12	12Cm	SDLC-1212-K-GS	SDLC-1212-KJ-GS	
12	15Cm	SDLC-1215-K-GS	SDLC-1215-KJ-GS	SDLC-1215-KPC-GS
12	17Cm	SDLC-1217-K-GS	SDLC-1217-KJ-GS	SDLC-1217-KPC-GS
12	20Cm	SDLC-1220-K-GS	SDLC-1220-KJ-GS	SDLC-1220-KPC-GS
12	25Cm	SDLC-1225-K-GS	SDLC-1225-KJ-GS	SDLC-1225-KPC-GS
14	15Cm	SDLC-1415-K-GS	SDLC-1415-KJ-GS	SDLC-1415-KPC-GS
14	17Cm	SDLC-1417-K-GS	SDLC-1417-KJ-GS	SDLC-1417-KPC-GS
14	20Cm	SDLC-1420-K-GS	SDLC-1420-KJ-GS	SDLC-1420-KPC-GS
14	25Cm	SDLC-1425-K-GS	SDLC-1425-KJ-GS	SDLC-1425-KPC-GS
14	30Cm	SDLC-1430-K-GS	SDLC-1430-KJ-GS	SDLC-1430-KPC-GS
14	35Cm	SDLC-1435-K-GS	SDLC-1435-KJ-GS	SDLC-1435-KPC-GS

Kit Contents

- (1) Step Tip Soft Dual Lumen Haemodialysis Subclavian Catheter
- Stiffener with Haemostatic Valve
- (1) J End Nitinol Guide Wire with Scaled Dispenser
- (2) Vessel Dilators
- (1) Guiding Syringe
- (1) Introducer Needle, Echogenic
- (1) Scalpel
- (2) Injection Cap
- (1) Unifix Catheter Tube Fixation Adhesive
- (2) Transparent Catheter Fixation Adhesive
- (1) Movable Fixation Wings (with Pre-curved catheter only)



1.2.4 Smart Dual Lumen Haemodialysis Catheter Kits

Ordering Information: STD Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
12	12 Cm	DLC-1212-KSM	DLC-1212-KJSM	DLC-1212-KPCSM
12	15 Cm	DLC-1215-KSM	DLC-1215-KJSM	DLC-1215-KPCSM
12	17 Cm	DLC-1217-KSM	DLC-1217-KJSM	DLC-1217-KPCSM
12	20 Cm	DLC-1220-KSM	DLC-1220-KJSM	DLC-1220-KPCSM
14	15 Cm	DLC-1415-KSM	DLC-1415-KJSM	DLC-1415-KPCSM
14	17 Cm	DLC-1417-KSM	DLC-1417-KJSM	DLC-1417-KPCSM
14	20 Cm	DLC-1420-KSM	DLC-1420-KJSM	DLC-1420-KPCSM

Kit Contents


- Dual lumen Smart Catheter
- J End Nitinol Guide Wire with Scaled Dispenser
- 2 Vessel Dilators
- Syringe
- Guiding Y Connector
- Introducer Needle, Echogenic
- Scalpel
- (2) Injection Cap
- Movable fixation Wing (with Pre-curved catheter only)

Ordering Information: GS Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
12	12 Cm	DLC-1212-KSM-GS	DLC-1212-KJSM-GS	DLC-1212-KPCSM-GS
12	15 Cm	DLC-1215-KSM-GS	DLC-1215-KJSM-GS	DLC-1215-KPCSM-GS
12	17 Cm	DLC-1217-KSM-GS	DLC-1217-KJSM-GS	DLC-1217-KPCSM-GS
12	20 Cm	DLC-1220-KSM-GS	DLC-1220-KJSM-GS	DLC-1220-KPCSM-GS
14	15 Cm	DLC-1415-KSM-GS	DLC-1415-KJSM-GS	DLC-1415-KPCSM-GS
14	17 Cm	DLC-1417-KSM-GS	DLC-1417-KJSM-GS	DLC-1417-KPCSM-GS
14	20 Cm	DLC-1420-KSM-GS	DLC-1420-KJSM-GS	DLC-1420-KPCSM-GS


Kit Contents

- Dual lumen Smart Catheter
- J End Nitinol Guide Wire with Scaled Dispenser
- 2 Vessel Dilators
- Guiding Syringe
- Introducer Needle, Echogenic
- Scalpel
- 2 Injection Cap
- Movable fixation Wing (with Pre-curved catheter only)



1.3. Triple Lumen catheters

1.3.1 Triple Lumen catheter kits



Ordering Information:
STD Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
12	15 Cm	TLC-1215-K	TLC-1215-KJ	TLC-1215-KPC
13	15 Cm	TLC-1315-K	TLC-1315-KJ	TLC-1315-KPC
13	17 Cm	TLC-1317-K	TLC-1317-KJ	TLC-1317- KPC
13	20 Cm	TLC-1320-K	TLC-1320-KJ	TLC-1320- KPC
13	24 Cm	TLC-1324-K	TLC-1324-KJ	TLC-1324-KPC
13	25 Cm	TLC-1325-K	TLC-1325-KJ	TLC-1325-KPC
14	15 Cm	TLC-1415-K	TLC-1415-KJ	TLC-1415- KPC
14	17 Cm	TLC-1417-K	TLC-1417-KJ	TLC-1417- KPC
14	20 Cm	TLC-1420-K	TLC-1420-KJ	TLC-1420- KPC

Kit Contents

- Triple Lumen Short Term Haemodialysis Catheter
- J End Nitinol Guide Wire with Scaled Dispenser
- 2 Vessel Dilators
- Syringe
- Guiding Y Connector
- Introducer Needle, Echogenic
- Scalpel
- 3 Injection Cap
- Movable fixation Wing (with Pre-curved catheter only)

Ordering Information:
STD Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
12	15 Cm	TLC-1215-K-GS	TLC-1215-KJ-GS	TLC-1215-KPC-GS
13	15 Cm	TLC-1315-K-GS	TLC-1315-KJ-GS	TLC-1315-KPC-GS
13	17 Cm	TLC-1317-K-GS	TLC-1317-KJ-GS	TLC-1317- KPC-GS
13	20 Cm	TLC-1320-K-GS	TLC-1320-KJ-GS	TLC-1320- KPC-GS
13	24 Cm	TLC-1324-K-GS	TLC-1324-KJ-GS	TLC-1324-KPC-GS
13	25 Cm	TLC-1325-K-GS	TLC-1325-KJ-GS	TLC-1325-KPC-GS
14	15 Cm	TLC-1415-K-GS	TLC-1415-KJ-GS	TLC-1415- KPC-GS
14	17 Cm	TLC-1417-K-GS	TLC-1417-KJ-GS	TLC-1417- KPC-GS
14	20 Cm	TLC-1420-K-GS	TLC-1420-KJ-GS	TLC-1420- KPC-GS

Kit Contents

- Triple Lumen Short Term Haemodialysis Catheter
- J End Nitinol Guide Wire with Scaled Dispenser
- 2 Vessel Dilators
- Guiding Syringe
- Introducer Needle, Echogenic
- Scalpel
- 3 Injection Cap
- Movable fixation Wing (with Pre-curved catheter only)

1.3.2. Smart Triple Lumen catheter kits

Ordering Information: STD Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
13	15 Cm	TLC-1315-KSM	TLC-1315-KJSM	TLC-1315-KPCSM
13	17 Cm	TLC-1317- KSM	TLC-1317- KJSM	TLC-1317- KPCSM
13	20 Cm	TLC-1320- KSM	TLC-1320- KJSM	TLC-1320- KPCSM
14	15 Cm	TLC-1415- KSM	TLC-1415- KJSM	TLC-1415- KPCSM
14	17 Cm	TLC-1417- KSM	TLC-1417- KJSM	TLC-1417- KPCSM
14	20 Cm	TLC-1420- KSM	TLC-1420- KJSM	TLC-1420- KPCSM

Kit Contents

- Triple Lumen Haemodialysis Subclavian Catheter
- J End Nitinol Guide Wire with Scaled Dispenser
- 2 Vessel Dilators
- Syringe
- Guiding Y Connector
- Introducer Needle, Echogenic
- Scalpel
- 3 Injection Cap
- Movable fixation Wing (with Pre-curved catheter only)

Ordering Information: GS Kit

Fr Size	Length	Subclavian	Jugular	Pre-curved
13	15 Cm	TLC-1315-KSM-GS	TLC-1315-KJSM-GS	TLC-1315-KPCSM-GS
13	17 Cm	TLC-1317-KSM-GS	TLC-1317-KJSM-GS	TLC-1317-KPCSM-GS
13	20 Cm	TLC-1320-KSM-GS	TLC-1320-KJSM-GS	TLC-1320-KPCSM-GS
14	15 Cm	TLC-1415-KSM-GS	TLC-1415-KJSM-GS	TLC-1415-KPCSM-GS
14	17 Cm	TLC-1417-KSM-GS	TLC-1417-KJSM-GS	TLC-1417-KPCSM-GS
14	20 Cm	TLC-1420-KSM-GS	TLC-1420-KJSM-GS	TLC-1420-KPCSM-GS

Kit Contents


- Triple Lumen Haemodialysis Subclavian Catheter
- J End Nitinol Guide Wire with Scaled Dispenser
- 2 Vessel Dilators
- Guiding Syringe
- Introducer Needle, Echogenic
- Scalpel
- 3 Injection Cap
- Movable fixation Wing (with Pre-curved catheter only)

2. Long Term haemodialysis Catheters

The Amecath Permthane Cuffed haemodialysis catheters are chronic, single or dual lumen, radiopaque, polyurethane catheters with a polyester cuff. The catheters are consisting of a catheter tube that have distally open, split or stepped tip and ends proximally with or without a connector hub. Catheters with hub are suitable for forward tunneling procedure while those without are suitable for retrograde tunneling. A female luer end is to be re-assembled to the catheter after the retrograde tunnelization.

Types and Differences between Amecath Permthane Dialysis Catheters:

Catheter type	Distal end tip design	Connect design	Proximal end design	Tunneling
Dual Lumen Catheter	stepped	with	Fixed only	Forward
Y Split Dual Lumen Catheter	split	with	Fixed only	Forward
Twin Cath Dual Lumen Catheter	stepped	without	re-assembly	retrograde
X Split Dual Lumen Catheter	split	without	re-assembly	retrograde
One Tube Catheter	rounded	without	re-assembly	retrograde



2.1. Dual Lumen long-term Haemodialysis catheter kit

Features and Properties

CATHETER FEATURES:

Tip design to guarantees soft contact and long patency.

PERMTHANE (Carbothane) material suitable for long use.

Carbothane provides excellent tensile strength.

Catheter have stepped tip and Available in straight and pre-curved shape. Extension lines stand repeated clamping and declamping.

Variety of sizes, lengths and shapes.

Ordering Information:

Fr. Size	Tip to Cuff Length (cm)	Total Length (cm)	Straight Codes	Pre-Curved Codes
12	16	21	PDLC-1216-K	PDLC-1216-KPC
12	19	24	PDLC-1219-K	PDLC-1219-KPC
14	19	24	PDLC-1419-K	PDLC-1419-KPC
14	23	28	PDLC-1423-K	PDLC-1423-KPC
14	27	32	PDLC-1427-K	PDLC-1427-KPC
14	31	36	PDLC-1431-K	PDLC-1431-KPC
14	35	40	PDLC-1435-K	PDLC-1435-KPC
14	50	55	PDLC-1450-K	PDLC-1450-KPC
16	19	24	PDLC-1619-K	PDLC-1619-KPC
16	23	28	PDLC-1623-K	PDLC-1623-KPC
16	27	32	PDLC-1627-K	PDLC-1627-KPC
16	31	36	PDLC-1631-K	PDLC-1631-KPC
16	35	40	PDLC-1635-K	PDLC-1635-KPC
16	50	55	PDLC-1650-K	PDLC-1650-KPC



Kit Contents

- (1) Permthane Long Term Dual Lumen Haemodialysis Catheter.
- (1) J End Nitinol Core Guide Wire
- (1) Stiffener with Cap
- (2) Vessel Dilators
- (1) 10cc Syringe
- (1) G18 x 7 cm Introducer Needle, Echogenic
- (1) Scalpel # 11
- (1) Valved Peelable Introducer Sheath (Peelable Introducer Sheath for 16Fr. Catheter)
- (1) Blunt Curved Tunneling Stylet
- (2) Injection Caps
- (2) Transparent Catheter Fixation Adhesive
- (1) Guiding Y Connector

2.2. Y-Split Dual Lumen long-term Haemodialysis catheter kit

Features and Properties

CATHETER FEATURES:

Tip design to guarantee soft contact and long patency.
 PERMTHANE (Carbothane) material suitable for long use.
 Carbothane provides excellent tensile strength.
 Catheter have split tip and Available in straight and pre-curved shape.
 Extension lines stand repeated clamping and declamping.
 Variety of sizes, lengths and shapes.

Ordering Information:

Fr. Size	Tip to Cuff Length (cm)	Total Length (cm)	Straight Codes	Pre-Curved Codes
12	16	21	PYDLC-1216-K	PYDLC-1216-KPC
12	19	24	PYDLC-1219-K	PYDLC-1219-KPC
14	19	24	PYDLC-1419-K	PYDLC-1419-KPC
14	23	28	PYDLC-1423-K	PYDLC-1423-KPC
14	27	32	PYDLC-1427-K	PYDLC-1427-KPC
14	31	36	PYDLC-1431-K	PYDLC-1431-KPC
14	35	40	PYDLC-1435-K	PYDLC-1435-KPC
14	50	55	PYDLC-1450-K	PYDLC-1450-KPC
16	19	24	PYDLC-1619-K	PYDLC-1619-KPC
16	23	28	PYDLC-1623-K	PYDLC-1623-KPC
16	27	32	PYDLC-1627-K	PYDLC-1627-KPC
16	31	36	PYDLC-1631-K	PYDLC-1631-KPC
16	35	40	PYDLC-1635-K	PYDLC-1635-KPC
16	50	55	PYDLC-1650-K	PYDLC-1650-KPC



Kit Contents

- (1) Permthane Long Term Y Split Dual Lumen Haemodialysis Catheter.
- (1) J End Nitinol Core Guide Wire
- (1) Stiffener with Cap
- (2) Vessel Dilators
- (1) 10cc Syringe
- (1) G18 x 7 cm Introducer Needle, Echogenic
- (1) Scalpel # 11
- (1) Valved Peelable Introducer Sheath (Peelable Introducer Sheath for 16Fr. Catheter)
- (1) Blunt Curved Tunnelling Stylet
- (2) Injection Caps
- (2) Transparent Catheter Fixation Adhesive
- (1) Guiding Y Connector

2.3. Twin lumen Haemodialysis Catheter “P2TC”

Features and Properties

The Twin lumen catheter design provides easy catheter insertion by a single tunneling procedure, through use of one outer wall configuration that ends with step tip. Catheter lumens are rounded in shape and this guarantees laminar high flow in each lumen.

2.3.1 P2TC with Integrated Female Luer Valve Kit

Ordering Information:

Fr. Size	Tip to Cuff Length (cm)	Total Length (cm)	Straight Codes
08	09	12	P2TC-0809-KV
08	12	15	P2TC-0812-KV
08	15	18	P2TC-0815-KV
10	09	12	P2TC-1009-KV
10	13	15	P2TC-1013-KV
10	15	18	P2TC-1015-KV
10	19	23	P2TC-1019-KV
12	24	28	P2TC-1224-KV
12	28	32	P2TC-1228-KV

Kit Contents

- (1) Permthane Long Term Pediatric Twin Lumen Catheter with Integrated Female Luer Valve
- (1) J End Nitinol Core Guide Wire
- (2) Vessel Dilators
- (1) 5 cc Syringe
- (1) G18 x 5 cm Introducer Needle, Echogenic
- (1) Scalpel #11
- (1) Movable Fixation Wing
- (1) Valved Peelable Introducer Sheath
- (1) Blunt Curved Tunneling Stylet
- (2) Repair Valves
- (2) Transparent Catheter Fixation Adhesive
- (1) Catheter Tube Fixation Adhesive (Unifix)
- (1) Guiding Y Connector

2.3.2 Twin lumen Haemodialysis Catheter Kit (Pediatric)

Ordering Information:

Fr. Size	Tip to Cuff Length (cm)	Total Length (cm)	Straight Codes	Pre-Curved Codes
08	10	13	P2TC-0810-K	P2TC-0810-KPC
08	12	15	P2TC-0812-K	P2TC-0812-KPC
10	12	15	P2TC-1012-K	P2TC-1012-KPC
10	15	18	P2TC-1015-K	P2TC-1015-KPC

Kit Contents

- (1) Permthane Long Term Pediatric Twin Lumen Catheter
- (1) J End Nitinol Core Guide Wire
- (2) Vessel Dilators
- (1) 5 cc Syringe
- (1) G18 x 5 cm Introducer Needle, Echogenic
- (1) Scalpel #11
- (1) Movable Fixation Wing
- (1) Valved Peelable Introducer Sheath
- (1) Blunt Curved Tunneling Stylet
- (2) Repair Valves
- (2) Transparent Catheter Fixation Adhesive
- (1) Catheter Tube Fixation Adhesive (Unifix)
- (1) Guiding Y Connector

2.3.3 Twin lumen Haemodialysis Catheter Kit (Adult)

Ordering Information:

Fr. Size	Tip to Cuff Length (cm)	Total Length (cm)	Straight Codes	Pre-Curved Codes
12	16	21	P2TC-1216-K	P2TC-1216-KPC
12	19	24	P2TC-1219-K	P2TC-1219-KPC
14	19	24	P2TC-1419-K	P2TC-1419-KPC
14	23	28	P2TC-1423-K	P2TC-1423-KPC
14	27	32	P2TC-1427-K	P2TC-1427-KPC
14	31	36	P2TC-1431-K	P2TC-1431-KPC
14	35	40	P2TC-1435-K	P2TC-1435-KPC
14	50	55	P2TC-1450-K	P2TC-1450-KPC
16	19	24	P2TC-1619-K	P2TC-1619-KPC
16	23	28	P2TC-1623-K	P2TC-1623-KPC
16	27	32	P2TC-1627-K	P2TC-1627-KPC
16	31	36	P2TC-1631-K	P2TC-1631-KPC
16	35	40	P2TC-1635-K	P2TC-1635-KPC
16	50	55	P2TC-1650-K	P2TC-1650-KPC
16	27	32	P2TC-1627-K	P2TC-1627-KPC
16	31	36	P2TC-1631-K	P2TC-1631-KPC
16	35	40	P2TC-1635-K	P2TC-1635-KPC
16	50	55	P2TC-1650-K	P2TC-1650-KPC



Kit Contents

- (1) Permthane Long Term Twin Lumen Catheter
- (1) J End Nitinol Core Guide Wire
- (2) Vessel Dilators
- (1) 10 cc Syringe
- (1) G18 x 7 cm Introducer Needle, Echogenic
- (1) Scalpel #11
- (1) Movable Fixation Wing
- (1) Valved Peelable Introducer Sheath (Peelable Introducer Sheath for 16Fr Catheter)
- (1) Blunt Curved Tunneling Stylet
- (2) Free Female Luer Line for Retrograde (Red for Arterial + Blue for Venous)
- (1) Flushing Connector
- (2) Slide Clamps
- (2) Transparent Catheter Fixation Adhesive
- (1) Catheter Tube Fixation Adhesive (Unifix)
- (1) Guiding Y Connector


2.4. X Split Dual Lumen Haemodialysis Catheter kit

Features and Properties

The split catheter design provides easy catheter insertion by a single tunneling procedure, through use of one outer wall configuration that ends with two fully independent catheter tubes which are joined together but can be separated by splitting. This enhances the free floating of each lumen.

Ordering Information:

Fr. Size	Tip to Cuff Length (cm)	Total Length (cm)	Straight Codes	Pre-Curved Codes
12	16	21	PXDLC-1216-K	PXDLC-1216-KPC
12	19	24	PXDLC-1219-K	PXDLC-1219-KPC
14	19	24	PXDLC-1419-K	PXDLC-1419-KPC
14	23	28	PXDLC-1423-K	PXDLC-1423-KPC
14	27	32	PXDLC-1427-K	PXDLC-1427-KPC
14	31	36	PXDLC-1431-K	PXDLC-1431-KPC
14	35	40	PXDLC-1435-K	PXDLC-1435-KPC
14	50	55	PXDLC-1450-K	PXDLC-1450-KPC
16	19	24	PXDLC-1619-K	PXDLC-1619-KPC
16	23	28	PYDLC-1623-K	PYDLC-1623-KPC
16	27	32	PXDLC-1627-K	PXDLC-1627-KPC
16	31	36	PXDLC-1631-K	PXDLC-1631-KPC
16	35	40	PXDLC-1635-K	PXDLC-1635-KPC
16	50	55	PXDLC-1650-K	PXDLC-1650-KPC



Kit Contents

- (1) Permthane Long Term X Split Dual Lumen Haemodialysis Catheter.
- (1) J End Nitinol Core Guide Wire
- (2) Vessel Dilators
- (1) 10cc Syringe
- (1) G18 x 7 cm Introducer Needle, Echogenic
- (1) Scalpel # 11
- (1) Movable Fixation Wing
- (1) Valved Peelable Introducer Sheath (Peelable Introducer Sheath for 16Fr. Catheter)
- (1) Blunt Curved Tunneling Stylet
- (2) Free Female Luer Line for Retrograde (Red for Arterial + Blue for Venous)
- (1) Flushing Connector
- (2) Slide Clamps
- (1) Catheter Tube Fixation Adhesive (Unifix)
- (2) Transparent Catheter Fixation Adhesive
- (1) Guiding Y Connector

3. Permanent Peritoneal Dialysis Catheter Kit

Features and Properties

Suitable for Percutaneous application.

Catheters for long term use, also available in two shapes, Straight and curled.

Both can be equipped with one or two cuffs.

Peritoneal catheters are available in many sizes, total lengths and tip to cuff lengths supplied with luer activated safety valve to ensure repeated hygienic exchanges.



3.1 Straight Permanent Peritoneal Dialysis Catheter Kit

Fr. Size	No. of Cuffs	Length (cm)	Codes
10	1	22	PSP-1022-1K
10	1	26	PSP-1026-1K
10	2	26	PSP-1026-2K
12	1	27	PSP-1227-1K
12	2	27	PSP-1227-2K
12	2	32	PSP-1232-2K
14	1	37	PSP-1437-1K
14	1	47	PSP-1447-1K
14	1	52	PSP-1452-1K
14	1	57	PSP-1457-1K
14	2	52	PSP-1452-2K
14	2	62	PSP-1462-2K
16	1	47	PSP-1647-1K
16	1	57	PSP-1657-1K
16	2	52	PSP-1652-2K
16	2	62	PSP-1662-2K

Kit Contents

- (1) Carbothane Peritoneal Catheter. Straight with one or two cuffs
- (1) J End Nitinol core 0.035" x 70 cm Guide wire
- (2) Vessel Dilators
- (1) 10 cc Syringe
- (1) G18 x 7 cm Introducer Needle, Echogenic
- (1) Scalpel # 11
- (1) Grip-Lock Peelable Introducer Sheath
- (1) Blunt Curved Tunneling Stylet 18cm.
- (1) Movable Fixation Wing
- (1) Haemostatic Valve
- (2) Transparent Catheter Fixation Adhesive
- (1) Guiding Y Connector
- (1) Free Female luer Line for Retrograde-Blue (For Cath. Size 10Fr only)
- (1) Female Luer (For Cath. Size 12, 14 and 16Fr)
- (1) Clamp (For Cath. Size 12, 14 and 16Fr)

3.2 Curled Permanent Peritoneal Dialysis Catheter Kit

Fr. Size	No. of Cuffs	Length (cm)	Codes
10	1	32	PCP-1032-1K
10	2	36	PCP-1036-2K
12	1	37	PCP-1237-1K
12	2	37	PCP-1237-2K
14	1	42	PCP-1442-1K
14	1	47	PCP-1447-1K
14	1	57	PCP-1457-1K
14	2	36	PCP-1436-2K
14	2	42	PCP-1442-2K
14	2	47	PCP-1447-2K
14	2	52	PCP-1452-2K
14	2	57	PCP-1457-2K
14	2	62	PCP-1462-2K
16	1	47	PCP-1647-1K
16	1	57	PCP-1657-1K
16	2	42	PCP-1642-2K
16	2	47	PCP-1647-2K
16	2	52	PCP-1652-2K
16	2	57	PCP-1657-2K
16	2	62	PCP-1662-2K

Kit Contents

- (1) Carbothane Peritoneal Catheter. Curled with one or two cuffs
- (1) J End Nitinol core 0.035" x 70 cm Guide wire
- (2) Vessel Dilators
- (1) 10 cc Syringe
- (1) G18 x 7 cm Introducer Needle, Echogenic
- (1) Scalpel # 11
- (1) Grip-Lock Peelable Introducer Sheath
- (1) Blunt Curved Tunneling Stylet 18cm.
- (1) Movable Fixation Wing
- (1) Haemostatic Valve
- (2) Transparent Catheter Fixation Adhesive
- (1) Guiding Y Connector
- (1) Free Female luer Line for Retrograde-Blue (For Cath. Size 10Fr only)
- (1) Female Luer (For Cath. Size 12, 14 and 16Fr)
- (1) Clamp (For Cath. Size 12, 14 and 16Fr)

- Established at 2001 for design, development and production of sterile dialysis catheter kits and accessories, Central venous kits and accessories, Urological catheters kits, stents, needles and accessories.
- 250 people are currently working at Ameco among them 5 doctors, 8 engineers, high caliber administrative, logistic and production .
- All manufacturing processes are under one roof. Extrusion, molding, over molding, tipping, punching, assembly, coating packaging and sterilization.
- Many Success stories for contact manufacturing and private labeling with major international companies.
- In 2002, Amecath received ISO13485 certificate in scope of design, development and production of catheter and kits, entailing that Amecath is applying full product quality assurance and Quality management system from design, production to service provision.
- In 2002, Amecath received EC mark (EC certificate) for all products, entailing that Amecath products can be freely placed in the market of 25 European countries and other countries with recognition of the European directives
- In 2005, Amecath received EC-Design examination for class III products (Central Venous Haemodialysis catheters for long term use), entailing that the device has passed the examination for safety and performance.
- In 2011, Amecath received Health Canada approval for central venous catheters, short term haemodialysis catheter and sheath introducer and urology double loop stents.
- In 2011, Amecath received Russian registration for central venous catheter and hemodialysis catheter, other products are under the registration process
- In 2012, Amecath received FDA approval for central venous catheter and short term hemodialysis catheter other products are currently under fda review.
- In 2013, Amecath received FDA approval for Dual lumen implantable long term haemodialysis catheters, sheath introducers and other products are currently under FDA review.
- In 2013, Amecath received EC-Design examination for class III products (Central Venous and peripherally inserted central catheters), entailing that the device has passed the examination for safety and performance
- In 2013, Amecath received health Canada Approval for long term implantable hemodialysis catheters

Manufactured By:

AMECO MEDICAL INDUSTRIES

Industrial Zone B4 Plot 119 east,
10th of Ramadan city - Egypt

Tel: +20.554.501321/ 2

Fax :+20.554.501224

Mail: support@amecathgroup.com

Website: www.amecathgroup.com

Authoized European Representative

Obelis s.a

Bd. Général Wahis 53
1030 Brussels, BELGIUM

Tel: +32.2. 732.59.54

Fax: +32.2.732.60.03

Email: mail@obelis.net

Website: www.obelis.net