

Cavidagel®

Sterile Wound Filler Hydrogel with Alginate



Description:

Cavidagel® is a Sterile Amorphous Hydrogel with alginate.

Mode of action:

- In case of a dry wound, **Cavidagel®** turns the dry wound into a moist one (hydrogel effect). The dry wound absorbs the water contained in **Cavidagel®** and promotes the dissolution of necrotic/damaged tissue.
- In case of a wet wound, the hydrocolloid particles of **Cavidagel®** are activated and consequently **Cavidagel®** absorbs excessive fluid from the wound and restores the water balance in the wound (hydrocolloid effect).
- Gel contains propylene glycol which acts as a humectant (moisturizer) which prevents the gel from drying out and improves its handling.

Features:

- Debriding action: **Cavidagel®** rehydrates necrotic tissue and aids debridement.
- Desloughing action: **Cavidagel®** absorbs slough and exudate without damaging fragile granulation tissue.
- Creates moist wound healing environment: **Cavidagel®** helps to create an optimal moist wound management environment which promotes healing, reduces the likelihood of scarring and helps soothe pain.
- Stable and counteract leakage and maceration: The viscous texture of the gel makes it stay in place, and even after absorption of debris and exudate the gel will remain cohesive and counteract leakage and maceration.
- Easy to remove: The gel is easy to remove due to high cohesion even after absorption of debris and exudate.


Intended Purpose:

- **Cavidagel®** is primarily intended for the treatment of necrotic and sloughy wounds, e.g. leg ulcers, pressure ulcers and non-infected diabetic foot ulcers.
- **Cavidagel®** may also be used on first and second degree burns.
- The gel may be used throughout the healing process to provide a moist wound healing environment in all types of wounds, except those mentioned under precautions.

Product	Code	Size	Units / Case	Case dimensions
Cavidagel Tube	CAVT15	15 gm	48	29x20x11cm
Cavidagel Tube	CAVT20	20 gm	48	
Cavidagel Tube	CAVT25	25 gm	48	
Cavidagel Tube	CAVT30	30 gm	48	24x20x11.5cm
Cavidagel Tube	CAVT50	50 gm	24	
Cavidagel Tube	CAVT75	75 gm	24	
Cavidagel Tube	CAVT85	85 gm	24	
Cavidagel Sachet	CAVS6	6 gm	1000	

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Cavidagel						
Issue date : 16-11-2016		Prepared by : Marina Medhat						
Colors : 		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	14/03/2019	Test changed by Rania Rafaat	Essential requirement					
3	13/02/2020	amending table of sizes	CRS request					
4	02-03-2020	Adding shipper case dimension	Marketing point of view					
5	16-11-2021	- Updating codes - Adding new size to table of sizes	File Maintenance					
6	03-08-2022	Updating product's photo	File Maintenance					