

Incifilm®



Description:

Incifilm®: is a thin transparent polyurethane film that covers and protects incision site during surgeries.

Features:

- **Transparent non-glare surface:** the incision site is clearly visible as the PU film is transparent without glare.
- **Tear resistant:** **Incifilm®** is thin but extremely strong, yet it can be cleanly cut with a scalpel.
- **Conformable and elastic:** **Incifilm®** is highly conformable and adapts well to body contours.
- **Highly breathable:** **Incifilm®** is permeable to water vapor, thus allows the skin to breathe and reduces risk of maceration.
- **Waterproof:** **Incifilm®** prevents water entry and reduces the risk of contamination.
- **Bacteria-proof and virus-proof:** **Incifilm®** provides a sterile working field as it acts as a barrier to bacteria and viruses.
- **Adheres right up to wound edge:** the drape provides good initial and permanent adhesion.
- **Easy to handle:** **Incifilm®** is provided with two colored liners that assist in the release of the drape onto the skin, which makes the application easier and saves time in patient preparation.

Intended Purpose:

Suitable for use as a mechanical barrier to facilitate all types of surgical operations.


Product code	Size	Pieces/Box	Boxes/Case
IN100100	10cmx10cm	10	4
IN100120	10cmx12cm	10	4
IN100140	10cmx14cm	10	4
IN100200	10cmx20cm	10	4
IN100250	10cmx25cm	10	4
IN100300	10cmx30cm	10	4
IN120250	12cmx25cm	10	4
IN150200	15cmx20cm	10	4
IN150250	15cmx25cm	10	4
IN150280	15cmx28cm	10	4
IN180180	18cmx18cm	10	4
IN200250	20cmx25cm	10	4
IN200300	20cmx30cm	10	4
IN200450	20cmx45cm	10	4
IN280300	28cmx30cm	10	4
IN280450	28cmx45cm	10	4
IN300300	30cmx30cm	10	4

InciFilm®

Product code	Size	Pieces/Box	Boxes/Case
IN300400	30cmx40cm	10	4
IN300450	30cmx45cm	10	4
IN300600	30cmx60cm	10	4
IN250350	25cmx35cm	10	4
IN350350	35cmx35cm	10	4
IN350700	35cmx70cm	10	4
IN350400	35cmx40cm	10	4
IN400400	40cmx40cm	10	4
IN400420	40cmx42cm	10	4
IN400450	40cmx45cm	10	4
IN400600	40cmx60cm	10	4
IN420450	42cmx45cm	10	4
IN250450	25cmx45cm	10	4
IN450450	45cmx45cm	10	4
IN450550	45cmx55cm	10	4
IN450600	45cmx60cm	10	4
IN450750	45cmx75cm	10	4
IN450900	45cmx90cm	10	4
IN500660	50cmx66cm	10	4
IN500800	50cmx80cm	10	4
IN560600	56cmx60cm	10	4
IN560800	56cmx80cm	10	4
IN560840	56cmx84cm	10	4
IN560850	56cmx85cm	10	4
IN350600	35cmx60cm	10	4
IN600600	60cmx60cm	10	4
IN600800	60cmx80cm	10	4
IN600900	60cmx90cm	10	4
IN600950	60cmx95cm	10	4
IN650750	65cmx75cm	10	4
IN800900	80cmx90cm	10	4
IN700900	70cmx90cm	10	4
IN850900	85cmx90cm	10	4
IN450650	45cmx65cm	10	4
IN150190	15cmx19cm	10	4
IN350500	35cmx50cm	10	4
IN600760	60cmx76cm	10	4
IN500760	50cmx76cm	10	4
IN250250	25cmx25cm	10	4
IN900900	90cmx90cm	10	4
IN250280	25cmx28cm	10	4
IN750750	75cmx75cm	10	4
IN830910	83cmx91cm	10	4
IN8501000	85cmx100cm	10	4
IN250380	25cmx38cm	10	4
IN440910	44cmx91cm	10	4
IN5001000	50cmx100cm	10	4
IN380410	38cmx41cm	10	4
IN200200	20cmx20cm	10	4
IN450820	45cmx82cm	10	4
IN820850	82cmx85cm	10	4
IN450500	45cmx50cm	10	4
IN400480	40cmx48cm	10	4
IN480500	48cmx50cm	10	4

- Amria - Alexandria - Egypt
 Tel. : 002 03 4500264 Fax : 002 03 4500263
 - Kafr-El-Zayat P.O.Box 38 Egypt
 Tel. : 002 040 2545643 Fax : 002 040 2545645
 E-mail : sales@pharmaplast-online.com

 **Pharmaplast**
 www.pharmaplast-online.com

Pharmaplast Amria 23512, Alexandria - Egypt www.pharmaplast-online.com		Brochure artwork template no: I-TM /AR/001-6 (Version 3)						
Customer : Pharmaplast		Product name : Incifilm						
Issue date : 15-10-2015		Prepared by : Marina Medhat						
Colors : 		Dimension : A4						
Issue Number	Revision Date	Description of Change	Rationale For Change	Amended by	Effective date	Approved by		
						Artwork manager	Marketing department	RA Department
1	N/A	New Issue	N/A					
2	22-04-2019	new text by Jostina Ramsis						
3	20-08-2019	adding new size in table of sizes	Marketing point of view					
4	25-05-2021	updating code	To be aligned with the new coding system					
5	05-03-2022	Adding new codes	file maintenance					
6	10-10-2022	Updating table	file maintenance					